

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property Organization  
International Bureau



(43) International Publication Date  
13 June 2002 (13.06.2002)


PCT

(10) International Publication Number  
WO 02/45598 A2

- (51) International Patent Classification<sup>7</sup>: A61B 17/22 (71) Applicant (for all designated States except US): LUMEND, INC. [US/US]; 400 Chesapeake Drive, Redwood City, CA 94063 (US).
- (21) International Application Number: PCT/US01/46871
- (22) International Filing Date: 5 December 2001 (05.12.2001) (72) Inventors; and (75) Inventors/Applicants (for US only): SPARKS, Kurt, D. [US/US]; 1618 Sand Hill Road, #406, Palo Alto, CA 94304 (US). EMERY, Jeffrey, L. [US/US]; 2104 Meadowview Place, San Mateo, CA 94401 (US). SEYBOLD, Brent, D. [US/US]; 2435 Armstrong Avenue, Santa Clara, CA 95050 (US). KUPIECKI, David, J. [US/US]; 3276 Market Street, San Francisco, CA 94114 (US). PINSON, C., Danielle [US/US]; 711 Eunice Avenue, Mountain View, California 94040 (US). MADSEN, Allen, W. [US/US]; 2878 Kilo Avenue, San Jose, CA 95124 (US). KELEHER, Michael, D. [US/US]; 34330 Dobson Way, Fremont, CA 94555 (US). SALINAS, Sergio [US/US]; 642 Scott Avenue, Redwood City, CA 94063 (US). CLARK, Benjamin, J. [US/US]; 3814 East Lake Way, Redwood City, CA 94062 (US). SELMON, Matthew, R. [US/US]; 99 Walnut Avenue, Atherton, CA 94027 (US).
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data:
- |            |                               |    |
|------------|-------------------------------|----|
| 60/251,756 | 5 December 2000 (05.12.2000)  | US |
| 60/255,729 | 14 December 2000 (14.12.2000) | US |
| 60/263,350 | 22 January 2001 (22.01.2001)  | US |
| 60/263,397 | 22 January 2001 (22.01.2001)  | US |
| 60/263,579 | 22 January 2001 (22.01.2001)  | US |
| 60/263,580 | 22 January 2001 (22.01.2001)  | US |
| 60/263,589 | 22 January 2001 (22.01.2001)  | US |
| 60/268,263 | 12 February 2001 (12.02.2001) | US |
| 60/301,537 | 27 June 2001 (27.06.2001)     | US |
| 60/329,936 | 17 October 2001 (17.10.2001)  | US |
- (74) Agents: GREGORY, Richard, L. et al.; Perkins Coie LLP, P.O. Box 2168, Menlo Park, CA 94026 (US).

[Continued on next page]

(54) Title: CATHETER SYSTEM FOR VASCULAR RE-ENTRY FROM A SUB-INTIMAL SPACE



(57) Abstract: A catheter system and corresponding methods are provided for accessing a blood vessel true lumen from a sub-intimal plane of the vessel. The catheter system includes visualization elements for determining the orientation of the true lumen with respect to the sub-intimal plane at an identified entry site from a position in the sub-intimal plane. The entry site is distal to a chronic total occlusion (CTO). The catheter system also includes a system for physically securing tissue of the sub-intimal plane at the entry site to the catheter system. The attaching system reduces or eliminates catheter float within the sub-intimal space. The catheter system further includes re-entry devices to establish and maintain a path from the sub-intimal plane back into the vessel true lumen.

WO 02/45598 A2

BEST AVAILABLE COPY



(81) **Designated States (national):** AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZM, ZW.

European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

**Published:**

— without international search report and to be republished upon receipt of that report

(84) **Designated States (regional):** ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM),

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

**CATHETER SYSTEM FOR VASCULAR RE-ENTRY FROM A**  
**SUB-INTIMAL SPACE**

**TECHNICAL FIELD**

5           The disclosed embodiments relate to catheter systems for crossing vascular occlusions.

**BACKGROUND**

          An interventional guide wire or other interventional device is often used in  
10   medical procedures that attempt to establish a pathway through a heavily stenosed or chronically occluded vessel. A chronically occluded vessel is referred to as containing a chronic total occlusion (CTO). During these procedures, the guide wire or device can only be of clinical benefit to establish vessel patency if it is advanced distally into the vessel true lumen.

15           At times during the process of advancing the guide wire or device through the stenosed vessel or CTO, and beyond the control of the operator, the guide wire or device may inadvertently enter into the wall of the vessel itself, i.e. the sub-intimal plane or space, or dissection plane. Once in this sub-intimal plane, it becomes difficult to navigate the guide wire or device through the sub-intimal  
20   tissue to re-gain access into the vessel true lumen at points distal to the occlusion. The layer of tissue that separates the vessel true lumen from the sub-intimal plane is typically in the range from 100 to 500 micrometers for vessels in the diameter range from 2mm to 4mm, and from 100 to 3000 microns, in the largest vessels of the body. The composition of the tissue may be such that no guide wire or  
25   interventional device currently on the market can re-access the true lumen.

**BRIEF DESCRIPTION OF THE FIGURES**

**Figure 1** is a catheter system of an embodiment including a rotational cutting element for re-entry and an imaging element for visualization.

**Figure 2** is a catheter system including a forward cutting element for re-entry, under an alternative embodiment, and an imaging element for visualization.

**Figure 3** is a catheter system including a reverse cutting element for re-entry, under yet another alternative embodiment, and an imaging element for visualization.

**Figure 4** is a catheter system including a lumen for receiving at least one of an optical fiber system, a radio frequency (RF) system, and a specialized guide wire, as another embodiment for re-entry, and an imaging element for visualization.

**Figure 5** is a catheter system embodiment including radio frequency (RF) electrodes or contacts for re-entry and an imaging element for visualization.

**Figure 6** is a catheter system embodiment including radio frequency (RF) electrodes or contacts of an alternative configuration for re-entry and an imaging element for visualization.

**Figure 7A** is a catheter system embodiment including separate radio frequency (RF) electrode ports and vacuum/guide wire ports, along with an imaging element for visualization.

**Figure 7B** is a catheter system embodiment including a lumen that accepts a fiber optic system or a guide wire, separate vacuum ports, and an imaging element for visualization.

**Figure 8** is a catheter system of an embodiment including tissue-holding skewers, a re-entry element or pierce tool, and an imaging element for visualization.

**Figure 9A** is a catheter system including a nosecone with an internal ramp to guide an internal cannula element or laser for re-entry, under an embodiment.

**Figure 9B** is a catheter system in an initial position prior to cannula deployment, under the embodiment of Figure 9A.

5        **Figure 9C** is a catheter system with a cannula deployed and a re-entry element advanced across sub-intimal tissue, under the embodiment of Figure 9A.

**Figure 9D** is a catheter system with a cannula advanced into a vessel true lumen, under the embodiment of Figure 9A.

**Figure 9E** is a catheter system following retraction of a re-entry element  
10        with the cannula maintained in the vessel true lumen, under the embodiment of Figure 9A.

**Figure 9F** is a catheter system with a guide wire advanced into a vessel true lumen, under the embodiment of Figure 9A.

**Figure 10A** is a catheter system including a nosecone with an internal ramp  
15        to guide a specialized guide wire for re-entry, under an alternative embodiment of Figures 9A and 9B.

**Figure 10B** is a catheter system including a nosecone, under the embodiment of Figure 10A, showing a specialized guide wire deploying from the internal ramp.

20        **Figure 10C** is a catheter system including a nosecone, under the embodiment of Figure 10A, showing a distal taper of the specialized guide wire deploying from the internal ramp.

**Figure 10D** is a catheter system including a nosecone, under the embodiment of Figure 10A, showing the distal taper of the specialized guide wire  
25        repositioning from the internal ramp through the nosecone slot.

**Figure 10E** is a catheter system including a nosecone, under the embodiment of Figure 10A, showing catheter removal from a treatment site following deployment of the specialized guide wire using the nosecone.

**Figure 11** is a catheter system including a nosecone with an internal ramp,  
5 under an embodiment, for use with a typical guide wire.

**Figure 12A** is a catheter system including a nosecone with an internal ramp and an internal slidably disposed tube in a retracted position, under an embodiment, for use with a typical guide wire, specialized guide wires, RF systems, and optical fiber systems.

10 **Figure 12B** is a catheter system including the nosecone with the internal ramp and internal slidably disposed tube, under the embodiment of Figure 12A, where the internal slidably disposed tube is in an extended position.

**Figure 13** is a catheter system including a curved distal catheter tip, under an embodiment, for guiding typical guide wires, specialized guide wires, optical  
15 fiber systems, and RF systems.

**Figure 14A** is a catheter system of an embodiment including a nosecone with an internal slidably disposed push ramp for guiding typical guide wires, specialized guide wires, optical fiber systems, and RF systems.

**Figure 14B** is a catheter system under the embodiment of Figure 14A  
20 showing the internal slidably disposed push ramp in a retracted position.

**Figure 15** is a catheter system including a dual lumen shaft that transitions distally to a single lumen shaft, under an embodiment.

**Figure 16** is a specialized guide wire including a distal end that is machined to provide cutting flutes, under an embodiment.

25 **Figure 17A** is an optical fiber system of an embodiment, capable of deployment in a catheter system, for delivering laser energy to a distal termination.

**Figure 17B** is an optical fiber system including a lumen, under an alternative embodiment of Figure 17A.

**Figure 18** is a slidably disposed element that translates within a catheter lumen and includes one or more electrodes capable of transmitting radio  
5 frequency energy, under an embodiment.

**Figure 19** is a rotational Intra-Vascular Ultrasound (IVUS) element including an integral distal tip for re-entry, under an embodiment.

**Figure 20** is a flow diagram for vascular re-entry from a sub-intimal space.

**Figure 21** is a flow diagram for identifying and determining orientation of a  
10 vessel true lumen using IVUS, under an embodiment.

**Figure 22** is a flow diagram for identifying and determining orientation of a vessel true lumen using OCT, under an embodiment.

**Figure 23** is a flow diagram for identifying and determining orientation of a vessel true lumen using optical fiber visualization systems, under an embodiment.

15 **Figure 24** is a flow diagram for identifying and determining orientation of a vessel true lumen using Doppler ultrasound systems, under an embodiment.

**Figure 25** is a flow diagram for identifying and determining orientation of a vessel true lumen using fluoroscopic marking systems, under an embodiment.

**Figure 26** is a flow diagram for securing sub-intimal tissue at a vessel re-  
20 entry site by evacuating fluid of the sub-intimal plane, under an embodiment.

**Figure 27** is a flow diagram for securing sub-intimal tissue at a vessel re-entry site using vacuum, under an embodiment.

**Figure 28** is a flow diagram for securing sub-intimal tissue at a vessel re-entry site using mechanical devices, under an embodiment.

25 **Figure 29** is a flow diagram for cutting through sub-intimal tissue into a true lumen, under an embodiment.

**Figure 30** is a flow diagram for piercing a pathway through sub-intimal tissue into a true lumen, under an embodiment.

**Figure 31** is a flow diagram for establishing a pathway through sub-intimal tissue into a true lumen, under an alternative embodiment.

5        **Figure 32** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under an embodiment.

**Figure 33** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a first alternative embodiment.

10       **Figure 34** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a second alternative embodiment.

**Figure 35** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a third alternative  
15       embodiment.

**Figure 36** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a fourth alternative embodiment.

**Figure 37** is a flow diagram for using a guide wire to establish a pathway  
20       through sub-intimal tissue into a vessel true lumen, under a fifth alternative embodiment.

**Figure 38** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a sixth alternative embodiment.



**Figure 39** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a seventh alternative embodiment.

**Figure 40** is a flow diagram for using a guide wire to establish a pathway  
5 through sub-intimal tissue into a vessel true lumen, under an eighth alternative embodiment.

**Figure 41** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a ninth alternative embodiment.

10 **Figure 42** is a flow diagram for establishing a path into a vessel true lumen using radio frequency (RF) energy, under an embodiment.

**Figure 43** is a flow diagram for establishing a path into a vessel true lumen using RF energy, under a first alternative embodiment.

**Figure 44** is a flow diagram for establishing a path into a vessel true lumen  
15 using RF energy, under a second alternative embodiment.

**Figure 45** is a flow diagram for establishing a path into a vessel true lumen using RF energy, under a third alternative embodiment.

**Figure 46** is a flow diagram for establishing a path into a vessel true lumen using laser energy, under an embodiment.

20 **Figure 47** is a flow diagram for establishing a path into a vessel true lumen using laser energy, under a first alternative embodiment.

**Figure 48** is a flow diagram for establishing a path into a vessel true lumen using laser energy, under a second alternative embodiment.

**Figure 49** is a flow diagram for establishing a path into a vessel true lumen  
25 using laser energy, under a third alternative embodiment.

**Figures 50A and 50B** are a flow diagram for establishing a path into a vessel true lumen using laser energy, under a fourth alternative embodiment.

**Figures 51A and 51B** are a flow diagram for establishing a path into a vessel true lumen using laser energy, under a fifth alternative embodiment.

5        **Figure 52** is a flow diagram for establishing a path into a vessel true lumen using laser energy, under a sixth alternative embodiment.

      In the drawings, the same reference numbers identify identical or substantially similar elements or acts. To easily identify the discussion of any particular element or act, the most significant digit or digits in a reference number  
10        refer to the Figure number in which that element is first introduced (e.g., element 902 is first introduced and discussed with respect to Figure 9).

      Figure numbers followed by the letters "A," "B," "C," etc. indicate either (1) that two or more Figures together form a complete Figure (e.g., Figures 50A and 50B together form a single, complete Figure 50), but are split between two or more  
15        Figures because of paper size restrictions, amount of viewable area within a computer screen window, etc., or (2) that two or more Figures represent alternative embodiments or methods under aspects of the invention.

      Unless described otherwise below, the construction and operation of the various blocks and components shown in the figures are of conventional design.  
20        As a result, such blocks need not be described in further detail herein, because they will be understood by those skilled in the relevant art. Such further detail is omitted for brevity and so as not to obscure the detailed description of the invention. Any modifications necessary to the blocks or components in the figures can be readily made by one skilled in the relevant art based on the detailed  
25        description provided herein.

The headings provided herein are for convenience only and do not necessarily affect the scope or meaning of the claimed invention.

#### DETAILED DESCRIPTION OF THE ILLUSTRATED EMBODIMENTS

5           A catheter-based system, or catheter system, is described for the purpose of gaining access to the true lumen of a blood vessel (coronary or peripheral artery or vein) from a space within the vessel wall itself, referred to herein as a sub-intimal plane, or dissection plane. Throughout this document, the various catheter embodiments are referred to as the re-entry catheter or catheter system.

10           The following description provides specific details for a thorough understanding of, and enabling description for, embodiments of the invention. However, one skilled in the art will understand that the invention may be practiced without these details. In other instances, well-known structures and functions have not been shown or described in detail to avoid unnecessarily obscuring the  
15 description of the embodiments of the invention.

          Unless described otherwise herein, the embodiments described herein are well known or described in detail in the above-noted and cross-referenced provisional patent applications. Indeed, much of the detailed description provided herein is explicitly disclosed in the provisional patent applications; most or all of the  
20 additional material of aspects of the invention will be recognized by those skilled in the relevant art as being inherent in the detailed description provided in such provisional patent applications, or well known to those skilled in the relevant art. Those skilled in the relevant art can implement aspects of the invention based on the detailed description provided in the provisional patent applications.

25           **Figures 1 through 19** show numerous different embodiments of catheter systems or platforms and associated components. For each system, the different

methods and steps described may be combined to construct unique embodiments.

Note that the various combinations of methods and steps yield additional embodiments of the catheter system. It is understood by those skilled in the art that these additional combinations/embodiments are intuitive in view of the

5 platforms presented herein. A description of the numerous embodiments now follows.

**Figure 1** is a catheter system 100 of an embodiment including a rotational cutting element 102 for re-entry and an imaging element 104 for visualization. In addition to the rotational or crown cutting element 102 and the imaging element  
10 104, the catheter system 100 includes an outer shaft 106 that houses at least one vacuum lumen 108 or port. The outer shaft 106 couples to a nosecone 110 that includes a guide wire lumen 112. A typical outside diameter of the outer shaft/nosecone is approximately 0.060 inches, while that of the cutting element is approximately 0.045 inches and that of the imaging element is approximately  
15 0.030 inches, but the embodiment is not so limited.

**Figure 2** is a catheter system 200 including a forward cutting element 202 for re-entry, under an alternative embodiment, and an imaging element 204 for visualization. In addition to the forward or beveled needle cutting element 202 and the imaging element 204, the catheter system 200 includes an outer shaft 206 that  
20 houses at least one vacuum lumen 208 or port. The outer shaft 206 couples to a nosecone 210 that includes a guide wire lumen 212.

**Figure 3** is a catheter system 300 including a reverse cutting element 302 for re-entry, under yet another alternative embodiment, and an imaging element 304 for visualization. In addition to the reverse cutting element 302 and the  
25 imaging element 304, the catheter system 300 includes an outer shaft 306 that

houses at least one vacuum lumen 308 or port. The outer shaft 306 couples to a nosecone 310 that includes a guide wire exit lumen 312.

Figure 4 is a catheter system 400 including at least one lumen 402 for receiving re-entry working elements 414 including optical fiber systems, rotational  
5 Intra-Vascular Ultrasound (IVUS) systems including those having a specialized distal tip, radio frequency (RF) systems, and specialized guide wires. The catheter system 400 includes a nosecone 410 coupled to a catheter outer body 406. The catheter outer body 406 includes a marker band for fluoroscopic visualization 422.

The nosecone 410 includes the lumen 402 for receiving any of a variety of  
10 working elements 414. The lumen 402 for receiving working elements terminates with an exit ramp 418 and a lateral exit port.

The nosecone 410 also includes a lumen 416 for receiving working elements including an imaging element 404 and a guide wire (not shown), for example. A visualization window 424 is included for use with the imaging element  
15 404. The lumens 402 and 416 may also serve as vacuum ports or pathways 420, but are not so limited. The nosecone 410 also includes a guide wire exit lumen 412 in a distal end.

Figure 5 is a catheter system embodiment 500 including radio frequency (RF) electrodes or contacts 502 for re-entry and an imaging element 504 for  
20 visualization. In addition to the RF electrodes 502 and the imaging element 504, the catheter system 500 includes an outer shaft 506 that houses at least one vacuum lumen 508 or port. The outer shaft 506 couples to a nosecone 510 that includes a distal guide wire exit lumen 512.

Figure 6 is a catheter system embodiment 600 including radio frequency  
25 (RF) electrodes or contacts 602 of an alternative configuration. In addition to the opposing RF electrodes 602 the catheter system 600 includes an imaging element

604 and at least one vacuum lumen 608 or port. The nosecone 610 includes a distal guide wire exit lumen 612, but other embodiments may not include the exit lumen 612.

**Figure 7A** is a catheter system embodiment 700 including a nosecone 710  
5 having separate ports 702 that accept radio frequency (RF) electrodes 712 for vessel tissue ablation. Further, the nosecone 710 includes a vacuum and guide wire port 704, or vacuum port. The vacuum port 704, besides delivering vacuum to the nosecone, accepts working elements including, for example, a guide wire 714. The nosecone also includes an imaging element 706 and a distal guide wire  
10 port 708.

**Figure 7B** is a catheter system embodiment 750 including a nosecone 760 with a lumen 754 that accepts working elements 764 including guide wires and fiber optic tissue systems, for example fiber optic tissue ablation systems. Further, the nosecone 760 includes separate vacuum lumens or ports 752 that deliver  
15 vacuum to the nosecone. The nosecone also includes an imaging element 756 and a distal guide wire port 758.

**Figure 8** is a catheter system 800 of an embodiment having a nosecone 810 that accommodates tissue-holding skewers 802, a re-entry element or pierce tool 804, and an imaging element 806 for visualization. The lumen 808 carrying  
20 the pierce tool 804 can be used to support introduction of a guide wire (not shown). The nosecone 810 also includes a guide wire exit lumen 812. The lumens 814 that accept the tissue holding skewers 802 also provide vacuum to the nosecone 810.

**Figure 9A** is a catheter system including a nosecone 902 with an internal  
25 ramp 904 to guide an internal cannula element and a re-entry laser system, under an embodiment. The nosecone 902 also includes a guide wire distal exit port 906,

and a cutout 908 to guide fluoroscopic alignment. A side view 912, top view 914 and front view 916 of the nosecone 902 are shown.

**Figure 9B** is a catheter system in an initial position prior to cannula deployment, under the embodiment of Figure 9A. The catheter shaft 930 and nosecone 902 is positioned within a sub-intimal plane 940 of the vessel wall 950. The cannula 920 is positioned in the catheter system at a position proximal to the exit ramp 904. The re-entry element 924 is advanced into the cannula 920. A variety of re-entry elements or devices 924 are deployable through the cannula 920, including typical guide wires, specialized guide wires (see **Figure 16** and the associated description herein), fiber optic systems (see **Figures 17A and 17B** and the associated description herein), RF electrode systems (see **Figure 18** and the associated description herein), and Intra-Vascular Ultrasound Systems (IVUS) (see **Figure 19** and the associated description herein). When the IVUS system is used for visualization, its visualization position 926 is shown. The catheter system of an embodiment also includes vacuum ports 922.

**Figure 9C** is a catheter system with the slidably disposed cannula 920 deployed from the nosecone 902 via the exit ramp 904. The re-entry element 924 advances across sub-intimal tissue 942, establishing a path into the vessel true lumen 960. **Figure 9D** shows the cannula 920 advanced into the vessel true lumen 960. **Figure 9E** is the catheter system following retraction of the re-entry element 924 with the cannula 920 maintained in the vessel true lumen 960, under the embodiment of Figure 9A. **Figure 9F** shows a guide wire 928 advanced into the vessel true lumen 960 through the cannula 920, following retraction of the re-entry element 924.

**Figure 10A** is a catheter system including a nosecone 1002 with an internal ramp to guide a specialized guide wire for re-entry, under an alternative

embodiment of Figures 9A and 9B. The nosecone 1002 of this embodiment has a modified guide wire exit port 1006 to accommodate the specialized guide wire or a specialized cannula. **Figures 10B-10E** show deployment of the specialized guide wire.

5           **Figure 10B** shows the specialized guide wire 1010 in an initial phase of deployment using the internal ramp 1012. The specialized guide wire 1010 of an embodiment includes at least one distal-taper section 1014. **Figure 10C** shows the specialized guide wire 1010 in a further state of deployment. **Figure 10D** shows that the nosecone slot 1016 allows the specialized guide wire 1010 to  
10   translate or reposition 1020 into the guide wire exit port 1006 as the distal taper section 1014 of the specialized guide wire 1010 reaches the top of the exit ramp 1012. **Figure 10E** shows that the nosecone 1002 and catheter system are retracted 1030 or removed from a treatment site following deployment of the specialized guide wire using the nosecone.

15           **Figure 11** is a catheter system 1100 including a nosecone 1102 coupled to a catheter shaft 1104. The nosecone 1102 includes an internal ramp 1106. A variety of working elements or devices 1108 are deployable using the internal ramp 1106, including typical guide wires, specialized guide wires (see **Figure 16** and the associated description herein), fiber optic systems (see **Figures 17A and**  
20   **17B** and the associated description herein), RF system components (see **Figure 18** and the associated description herein), and IVUS (see **Figure 19**). The catheter system 1100 further includes a vacuum lumen 1110 and at least one region 1112 housing fluoroscopic visualization elements.

**Figure 12A** is a catheter system 1200 including a nosecone 1202 with an  
25   internal ramp 1204 and an internal slidably disposed tube 1206, under an alternative embodiment of **Figures 9A, 9B, and 10A-10E**. The position of the



internal tube 1206 is controllable to aid in steering a working element 1208 deployed via the catheter system 1200. When the internal tube 1206 is in this retracted position, the working element 1208 is deployed at a shallower deployment angle relative to a longitudinal axis of the catheter system 1200.

5 **Figure 12B** shows the catheter system 1200 when the internal tube 1206 is in an extended position. Extension of the internal tube 1206 results in deployment of the working element 1208 at deployment angles that are progressively more acute. A variety of working elements or devices 1208 are deployable using the internal ramp 1204, including typical guide wires, specialized guide wires (see **Figure 16**  
10 and the associated description herein), fiber optic systems (see **Figures 17A and 17B** and the associated description herein), RF system components (see **Figure 18** and the associated description herein), and IVUS (see **Figure 19** and the associated description herein).

**Figure 13** is a catheter system 1300 including a curved distal catheter tip  
15 1302, under an embodiment. The curved distal tip 1302 is coupled to a distal coil 1304 that, in one embodiment, is formed from platinum. The distal coil 1304 may include a visualization window 1306, but is not so limited. The distal coil 1304 couples to a braided catheter shaft 1308, but may be used with various types of catheter shafts known in the art. The catheter system 1300 can be used for  
20 guiding various working elements including typical guide wires, specialized guide wires, optical fiber systems (see **Figures 17A and 17B** and the associated description herein), and IVUS (see **Figure 19** and the associated description herein).

**Figures 14A and 14B** show a catheter system 1400 including a catheter  
25 shaft 1401 having a nosecone 1402 with an internal slidably disposed push ramp 1404 coupled to an internal push tube 1406. Extension of the internal push ramp

1404 (**Figure 14A**) helps in directing a working element 1410 to a re-entry site.

The nosecone 1402 further includes a nosecone slot 1408 that allows the working element 1410 to translate or reposition into the guide wire exit port 1412 as the internal ramp 1404 is retracted (**Figure 14B**). A variety of working elements or devices 1410 are deployable using the internal push ramp 1404, including typical guide wires, specialized guide wires (see **Figure 16** and the associated description herein), fiber optic systems (see **Figures 17A and 17B** and the associated description herein), RF system components (see **Figure 18** and the associated description herein), and IVUS (see **Figure 19** and the associated description herein).

**Figure 15** is a catheter system 1500 including a catheter shaft having dual lumens 1502 and 1504 in a proximal region, under an alternative embodiment of **Figure 13**. In a distal region, the two lumens 1502 and 1504 merge to form a single lumen 1506. A variety of working elements or devices are deployable using this catheter system 1500, including typical guide wires, specialized guide wires (see **Figure 16** and the associated description herein), fiber optic systems (see **Figures 17A and 17B** and the associated description herein), RF system components (see **Figure 18** and the associated description herein), and IVUS (see **Figure 19** and the associated description herein).

**Figure 16** is a specialized guide wire 1600, under an embodiment. The guide wire 1600 includes a distal end that is machined to provide cutting flutes 1602, under an embodiment. In an alternative embodiment, the distal tip may be processed with a generally abrasive surface.

**Figure 17A** is an optical fiber system 1700 of an embodiment. This optical fiber system 1700 can be deployed, using any of the catheter systems described herein, to deliver laser energy to a distal termination. The optical fiber system

1700 includes a fiber optic core 1702 surrounded by an outer sheath 1704. The sheath 1704 can be formed, for example, from polyimide or polyethylene, but is not so limited. The fiber optic system 1700 terminates at the distal end in either a normal configuration 1706 or an angled configuration 1708, relative to the longitudinal axis, but any mode of termination known in the art may be used. The optical fiber system 1700 of one embodiment has an outside diameter of approximately 0.010 to 0.018 inches, but the diameter can vary depending on the planned application.

Figure 17B is an optical fiber system 1750 including a lumen 1752, under an alternative embodiment of Figure 17A. The lumen 1752 is surrounded by a polymer encasement 1754. The encasement 1754 of one embodiment is formed, for example, from nylon or polyethylene, but is not so limited. The encasement 1754 includes individual optical fibers 1756, where the number of optical fibers 1756 varies in accordance with planned applications. The individual fibers can terminate in any mode known in the art. The optical fiber system 1750 of one embodiment has an outside diameter of approximately 0.020 to 0.030 inches, but the diameter can vary with particular applications.

Figure 18 is a radio frequency (RF) element 1800 including one or more electrodes 1802 capable of transmitting RF energy, under an embodiment. The RF electrodes 1802 are housed in a polymer sheathing 1804 that translates within a catheter lumen (not shown). The polymer sheathing 1804 is formed from polyethylene or nylon, but other materials may be used. The RF element 1800 can be used with the catheter system embodiments described herein.

Figure 19 is a rotational IVUS element 1900 including an integral distal tip 1902 for re-entry, under an embodiment. The distal tip 1902 is coupled to a torqueable drive train 1904. The element 1900 further includes an imaging

element 1906 and an imaging transducer 1908. The element 1900 is encased in a sleeve 1910 that, in one embodiment, is formed of polyethylene. The outside diameter of the element is approximately 0.030 inches, but is not so limited.

All catheter systems presented herein are delivered to the vascular site via tracking over a conventional guide wire. In some instances the guide wire is removed and other catheter elements are advanced during the course of a procedure involving the catheter system while, or in other instances, the guide wire remains in the catheter throughout the procedure.

Figure 20 is a flow diagram for vascular re-entry from a sub-intimal space.

In general, the process of re-entry from a sub-intimal plane into a vessel true lumen of an embodiment is described herein using three steps, with numerous methods and embodiments described under each step. Fundamentally, any methods from any steps can be combined to formulate a valid sequence and basis for a catheter embodiment to describe the overall procedures. The three steps include:

Step 1: Identify and determine the orientation of the vessel true lumen with respect to the sub-intimal plane. Approaches are presented under this step including the use of catheter system on-board guidance, and external guidance. Further, five methods of visualization are described.

Step 2: Physically secure the sub-intimal tissue at the re-entry site, to enable a method of re-entry into the true lumen, as described in Step 3. Three methods are described herein for securing the sub-intimal tissue.

Step 3: Establish a re-entry path from the sub-intimal plane into the vessel true lumen. Six methods of vessel re-entry are described below.

The steps are now described in detail, including the associated methods and embodiments.

Step 1: Identify and determine the orientation of the vessel true lumen with respect to the sub-intimal plane

As the catheter is positioned within a sub-intimal plane, the re-entry mechanism is orientated towards the vessel true lumen. When the catheter is properly aligned, the re-entry mechanism directly faces the sub-intimal tissue that separates the dissection plane from the vessel true lumen.

Method 1 under Step 1: Intra-Vascular Ultrasound (IVUS) (Figures 1-15)

Figure 21 is a flow diagram for identifying and determining orientation of a vessel true lumen using IVUS, under an embodiment. Two IVUS systems are readily available as stand-alone devices that may serve as an element within the Re-Entry Catheter. One system, manufactured by Boston Scientific Corporation, utilizes an ultrasound element, or crystal, which is mounted at the distal end of a rotational shaft. This shaft is rotated at a specified speed while the crystal is excited by electrical signals. The crystal produces acoustic wavefronts, and the reflection of the acoustic wavefronts by tissue types of varying density are received by the crystal. Algorithms decipher the reflected acoustic signals, and the system provides a cross sectional view of the tissue surrounding the sub-intimal plane. This method easily resolves surrounding tissue and aids in identifying the vessel true lumen.

Another IVUS system, manufactured by Endosonics utilizes an array of individually mounted crystals or elements in a fixed circumferential orientation at the distal end of a small catheter shaft. This system does not rotate. Rather, each

crystal is sequentially excited by an electrical signal and each crystal also receives the reflected acoustic signal. Algorithms decipher the reflected signals, and the system constructs a cross sectional image of the tissue surrounding the crystal network. This method resolves surrounding tissue and helps identify the vessel true lumen.

Typical outer diameter dimensions for IVUS systems are on the order of approximately 0.030 inches, but are not so limited.

Method 2 under Step 1: Optical Coherence Tomography (OCT) (Figures 1-

15)

Figure 22 is a flow diagram for identifying and determining orientation of a vessel true lumen using OCT, under an embodiment. In general, the OCT system delivers infrared (IR) light into tissue from the distal end of a rotating optical fiber. Delivery of the light into the tissue is accomplished by terminating the optical fiber at an angle to achieve internal reflection of the light at an approximate right angle to the central axis of the optical fiber. The fiber also receives reflected light from various tissue types. While the reflected OCT signals comprise light waves rather than acoustic waves, the reflected signals are deciphered and a cross sectional image is produced of the tissue surrounding the distal tip of the optical fiber, similar to the approach described in Method 1 under Step 1 above. While this technology has yet to be released to the market as a tool to be used within the vascular space, OCT is a proven technology that provides image resolution 5 to 25 times greater than current ultrasound embodiments.

Further, OCT technology is not able to produce an image through blood because it uses infrared light. However, the system of an embodiment evacuates fluids and/or blood from the sub-intimal plane in an embodiment, as described

below. Thus, with the removal of fluids and/or blood from the sub-intimal plane, the utilization of infrared light wavelengths allows the OCT technology to create an image of the surrounding tissue and identify the true lumen. Further, the wavelength of light may be modified so as to be able to pass through blood and still produce an image of the surrounding tissue. Typical outer diameter dimensions for the OCT system are on the order of approximately 0.015 inches, but are not so limited.

**Method 3 under Step 1: Optical Fiber Visualization System (Figures 1- 15)**

10       **Figure 23** is a flow diagram for identifying and determining orientation of a vessel true lumen using optical fiber visualization systems, under an embodiment. This method uses an optical fiber system for visualization as in Method 2 under Step 1 above, except that the optical fiber system does not continuously rotate. A broad-spectrum light is applied to the proximal end of the optical fiber. The light is reflected at the distal terminal end of the fiber into surrounding body tissue by terminating the optical fiber at an approximate right angle. Specific wavelengths of light are absorbed and/or transmitted and/or reflected based upon the type of tissue exposed to the light. The reflected light signals are processed to produce a spectral graph. This method, therefore, incrementally rotates the fiber, sending and receiving broadband light signals, until the specific absorption peaks are received that signal the presence of the true lumen. Typical outer diameter dimensions of the fiber optic system may range from 0.010 to 0.020 inches, but are not so limited.

25       This method can be implemented using three techniques that are now described.

Technique 1 (Method 3 under Step 1):

The absorption peak of hemoglobin is sought, which indicates the presence of a blood pool, i.e. the true lumen. Note that in the case where re-entry of an artery is desired, the absorption peak sought is that of oxygenated hemoglobin, which is a double peak signal having peaks at known wavelengths in the spectrum. Since large veins also are in close proximity to arteries, it would be important to distinguish between the single absorption peak of de-oxygenated hemoglobin, and the desired double peak of oxygenated hemoglobin. Once the double absorption peaks of de-oxygenated are identified, the direction of the true lumen is determined.

Technique 2 (Method 3 under Step 1):

Local injection into the vessel true lumen with an agent having a characteristic absorption peak can selectively tag blood cells only in the vessel true lumen. Various tagging agents approved by the Food and Drug Administration (FDA) can be used to tag blood cells. Tagging the blood cells in the vessel true lumen is accomplished from a distal entry point to the vessel.

As an example, consider an occlusion in the right coronary artery of the coronary vasculature. The re-entry catheter described herein is advanced into the dissection plane. By definition, the vessel true lumen at this time in the procedure has not been accessed. However, many times collateral vessels from the left coronary vasculature will branch to the right coronary artery at various locations, and often at connection points which are distal to typical occlusion locations in the right coronary artery. Thus, by injecting this agent into the left coronary vasculature, some of the agent will travel through the collateral vessels and feed into the right coronary artery, distal to the occlusion, i.e. the true lumen of the



vessel the re-entry catheter is attempting to access. This procedure uses efficient timing of the injection of the agent and interpretation of the absorption peaks, because the agent will clear from the distal portion of the right coronary artery typically within 5-10 seconds. Note that this example is easily applicable to a  
5 blockage in the left arterial system, with collateral connections from the right coronary artery.

Technique 3 (Method 3 under Step 1):

A systemic injection with an agent that tags red blood cells is used. In this  
10 case, the agent tags all viable red blood cells in the entire vascular system, both arterial and venous. The agent is chosen to provide an identifiable absorption peak. One advantage of this method is that a continuous generation of spectral data can be performed since the agent bond to the blood cells has a long half-life.

15 Method 4 under Step 1: Doppler Ultrasound (Figures 1-8, 15)

Figure 24 is a flow diagram for identifying and determining orientation of a vessel true lumen using Doppler ultrasound systems, under an embodiment. Doppler ultrasound is also used to identify blood flow in the vessel true lumen, under an embodiment. The fundamental basis of this technique is not unlike that  
20 used to generate weather radar maps. The Doppler method emits acoustic energy from a transducer/receiver mounted in the catheter. The transducer/receiver subsequently recognizes and measures phase shifts in these acoustic signals that are reflected off of moving, formed blood components, e.g. red blood cells and/or white blood cells in the vessel true lumen, thereby establishing an orientation.  
25 Details of the Doppler method are understood by those knowledgeable in this art.

Typical outer diameter dimensions of the doppler system are approximately 0.020 inches.

Note that the visualization hardware referenced in Methods 1 through 4 under Step 1 may be an integral component of the embodiments of Figures 7, 8, and 15, i.e., not removed from the catheter during the procedure. However, in the embodiments of Figures 1 through 6 and 9 through 14 this hardware is used as part of the catheter system at the beginning of the procedure to determine orientation of the vessel true lumen, and then removed to allow the introduction of other elements used to conduct the procedure.

10

Method 5 under Step 1: Fluoroscopic Marking System (Figures 1- 15)

Figure 25 is a flow diagram for identifying and determining orientation of a vessel true lumen using fluoroscopic marking systems, under an embodiment. This method describes the use of a fluoroscopic marking system at the distal end of the catheter to identify the location of the re-entry mechanism, with respect to the vessel true lumen. One implementation of this method takes advantage of the catheter nosecone itself which may be fabricated from a fluoroscopic material, e.g. stainless steel, platinum, or gold coated ceramic. For example the side port of the nosecone, which identifies the re-entry direction of the catheter, should be readily visible under fluoroscopy. Alternatively, a marking cutout may be placed within the proximal section of the nosecone, the design of which would indicate the position of the re-entry mechanism of the catheter.

Step 2: Physically secure the sub-intimal tissue at the re-entry site

In the formation of a sub-intimal plane, extra-vascular fluid as well as blood may collect in this space, and the sub-intimal space may grow in volume. A re-

entry catheter advanced into the sub-intimal space can thus "float" within this space. This can result in the catheter's inability to establish a "purchase" on the sub-intimal tissue that separates the catheter from the vessel true lumen. Thus, it is procedurally important to evacuate this volume of fluid so as to allow the re-entry portion of the catheter to be placed in direct, intimate contact with the sub-intimal tissue. This increases the likelihood of the re-entry mechanism successfully establishing a re-entry path to the vessel true lumen.

10      Method 1 under Step 2: Evacuating Fluid of Sub-Intimal Plane (Figures 1-15)

Figure 26 is a flow diagram for securing sub-intimal tissue at a vessel re-entry site by evacuating fluid of the sub-intimal plane, under an embodiment. This method describes the action of evacuating the volume of fluid that is typically contained within the sub-intimal plane. The action of evacuating the sub-intimal plane of all extra-vascular fluid and blood, can subsequently "lock" the sub-intimal tissue onto the surface of the catheter. Thus, the sub-intimal tissue that separates the sub-intimal plane from the true lumen is held in position on the surface of the re-entry catheter, facilitating a chosen method of creating a re-entry pathway back to the true lumen.

20      More specifically, the embodiments described herein and shown in the Figures include a distal nosecone having a sideport or catheter distal termination with ports which may be used to translate a vacuum to within the sub-intimal space. Note that in some embodiments the distal guide wire port serves to communicate vacuum to within the sub-intimal space. Further, additional  
25      evacuation/vacuum ports of various sizes and shapes may be positioned along the distal portion of the catheter which would reside in the sub-intimal plane.

**Method 2 under Step 2: Application of Vacuum (Figures 1 - 3, 5 - 12, 14)**

**Figure 27** is a flow diagram for securing sub-intimal tissue at a vessel re-entry site using vacuum, under an embodiment. This method describes the further application of vacuum such that the sub-intimal tissue that separates the sub-intimal plane from the vessel true lumen is invaginated within a distal portion of the catheter. Once this tissue is contained within the interior of the catheter itself, various methods may be subsequently employed to establish a physical pathway through it and back to the vessel true lumen, as detailed below. One advantage of this method is that all re-entry techniques may be performed within the catheter, reducing risk to any surrounding vascular tissue and inadvertent perforation of the vessel wall to the pericardial space.

**Method 3 under Step 2: Mechanical Means (Figure 8)**

**Figure 28** is a flow diagram for securing sub-intimal tissue at a vessel re-entry site using mechanical devices, under an embodiment. This method describes mechanically securing the sub-intimal tissue prior to employment of a re-entry mechanism to gain access into the vessel true lumen, per Step 3 below. Physical or mechanical securing of the sub-intimal tissue may be accomplished by a variety of embodiments. Mechanical means including tweezing or forcep action, pinch rollers, or skewers are used to grab and secure the sub-intimal tissue. The tweezing- or forcep-type mechanical system can be used to simply hold the tissue secure against the nosecone, or grab and pull the sub-intimal tissue within the nosecone. Skewers may be used to pierce the sub-intimal tissue and hold it intimately against the catheter, or grab and pull the sub-intimal tissue within the nosecone. The skewer tip could remain embedded in the sub-intimal tissue, or could advance through the tissue such that the skewer tip pierces in and out of the

tissue and the skewer is directed into a receiving port at the opposite section of the nosecone. In this last configuration, the tissue is held captive on the skewer.

Step 3: Establish a re-entry path from the sub-intimal plane into the vessel

5 true lumen

This section describes the methods to establish the re-entry pathway. As each method may apply to various ones of the embodiments herein, numerous combinations of method/catheter platform are presented.

10 Method 1 under Step 3: Cutting Sub-Intimal Tissue (Figures 1- 3)

This method describes the cutting of a pathway through the sub-intimal tissue into the true lumen. Three embodiments for this method include a catheter shaft with a nosecone termination which houses the visualization element (Step 1, Methods 1-4) and the cutting member as shown in the referenced figures, and  
15 further described below. The nosecone has a sideport opening which allows for visualization of the vascular area (per Step 1, Methods 1 through 4) and identification of the orientation of the vessel true lumen. Fluoroscopic features of the nosecone can also be used in conjunction with the primary visualization method to facilitate alignment to the true lumen. The nosecone also has a distal  
20 end port that allows the catheter to be tracked in a co-linear fashion over a conventional guide wire to the chosen vascular site.

Contained within the outer shaft is an internal shaft to which a distal cutting element is attached. The catheter outer shaft and the internal shaft may be fabricated using any number of methods know in the art. An example would be  
25 polymer lamination onto a stainless steel braided tube. Alternatively, either shaft may be fabricated as described in United States Patent Application Number

09/984,498, filed October 16, 2001. The cutting member may be fabricated starting with a stainless steel hypotube, and forming the appropriate cutting features, e.g. serrations (Figure 1), needle point (Figures 2 and 3), or a sharpened conical termination (not shown). The cutting features may be formed using  
5 standard machining methods or electronic discharge machining (EDM).

Representative dimensions of a cutting element may range from approximately 0.030 to 0.040 inches in diameter. The nosecone may be fabricated of similar materials, using similar fabrication methods.

Figure 29 is a flow diagram for cutting through sub-intimal tissue into a true  
10 lumen, under an embodiment. Procedurally, prior to introducing the catheter into the vasculature, the imaging element is removed from the catheter. The inner shaft/cutting element is first positioned such that it covers the nosecone sideport. For the embodiments of Figures 1 and 2, the cutting element is advanced fully distal such that the distal termination of the cutting element is "garaged" within the  
15 distal receiving section of the nosecone. In the embodiment of Figure 3, the cutting element is rotated until it is positioned opposite the nosecone sideport, e.g., the sideport is covered by the wall of the hypotube opposite the cutting feature. In these configurations, the cutting element covers the nosecone sideport and presents an uninterrupted surface that can be effectively tracked through the  
20 vasculature. The catheter may then be introduced onto a guide wire and tracked to the appropriate vascular location.

Once the vascular site has been reached, the guide wire is removed and the imaging element is introduced and advanced until it is positioned at the nosecone sideport. The cutting element is retracted, exposing the sideport. The  
25 imaging element is activated and the sideport is rotationally positioned towards the vessel true lumen. Vacuum is applied to the interior of the catheter outer shaft per

Step 2, Method 2, thus beginning the process to invaginate the sub-intimal tissue into the nosecone sideport. While confirming that the sideport remains directed towards the vessel true lumen, the imaging element is retracted, allowing more space for the sub-intimal tissue to further invaginate into the nosecone sideport.

5 Note that under this configuration, visualization of the vessel true lumen is no longer possible since the imaging element is retracted proximal to the sideport.

While maintaining vacuum, the cutter is then advanced, distally for the embodiments of Figures 1 and 2, and proximally for the embodiment of Figure 3. In the embodiment of Figure 1, the cutting element may also be rotated to facilitate

10 the cutting action. In this process, the actuation of the cutting element mechanically traps and compresses the tissue within the sideport, allowing the cutting features to propagate an incision through the sub-intimal tissue, and through to the vessel true lumen.

Note that an alternative to the process described above is to not retract the

15 visualization element prior to the cutting action. Thus the visualization element is positioned at the nosecone sideport, and cutting is performed while simultaneously viewing the alignment to the vessel true lumen. Also note that in this alternative technique, the presence of the visualization element in the nosecone prevents maximum invagination of tissue into the nosecone. Thus the first technique allows

20 cutting through thicker sub-intimal tissue which separates the sub-intimal plane from the vessel true lumen.

Upon establishment of a pathway into the vessel true lumen, the vacuum still applied at the nosecone will aspirate blood from the vessel true lumen, through the pathway in the sub-intimal tissue and into the catheter shaft, ultimately

25 reaching the proximal end of the catheter. The vacuum will also be lost. These

two events would indicate that a pathway has been successfully established into the vessel true lumen.

Once a pathway has been established into the vessel true lumen, the imaging element is removed, the cutting element is once again positioned to cover the nosecone sideport, and a guide wire is advanced through the catheter to exit the nosecone endport. The guide wire is thus manipulated through the pathway into the vessel and the catheter is removed from the vasculature.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.060 to 0.070 inches; cutting element outside diameter is approximately 0.040 to 0.050 inches; and imaging element outside diameter is as described above.

#### Method 2 under Step 3: Piercing (Figures 8 and 9)

This method describes the general piercing of a pathway through the sub-intimal tissue into the true lumen. This can be approached under a first embodiment where the sub-intimal tissue is held onto the surface of the catheter via vacuum (Figure 9), or under a second embodiment where the sub-intimal tissue is invaginated into the catheter via vacuum (Figure 8). These two embodiments are now described.

#### Embodiment 1 (Method 2 under Step 3) (Figure 9)

This embodiment includes a nosecone or molded catheter termination attached to the distal end of the catheter, and an internal slidably disposed actuating cannula.

The catheter shaft may be any of a number of catheter shafts known in the art. The nosecone includes a side exit port and a distal end port coupled via a slot



which, and as will be described below, allows the guide wire to move from the side port into the distal end port when the catheter is retracted proximally over the guide wire. The cannula is guided out of the nosecone sideport via two internal exit ramps, one on either side of the cannula. The distal end port of the nosecone  
5 allows the catheter to be tracked in a co-linear fashion over a standard coronary guide wire.

Figure 30 is a flow diagram for piercing a pathway through sub-intimal tissue into a vessel true lumen under an embodiment. Procedurally, a guide wire is placed in the sub-intimal space of the target vasculature such that the guide wire  
10 distal end is located distal to the occluded area of the vessel. The cannula is retracted to a position proximal to the exit ramp so that the cannula exit port is co-linear with the inner diameter of the nosecone. This configuration allows the proximal end of the guide wire to be passed through the nosecone distal end port. The cannula, the catheter shaft, and thus the catheter may be tracked over the  
15 guide wire to the vascular site.

The catheter is then aligned to the vessel true lumen. This is accomplished per Step 1, Methods 1 through 4, or Step 1, Method 5. Figure 9 illustrates both configurations. In the case using Step 1, Methods 1 through 4, the guide wire is retracted from the catheter, and the visualization element is advance into the  
20 nosecone. The element is activated at the nosecone side port and the side port is rotated to face the vessel true lumen. The visualization element is removed, and the cannula and guide wire are re-introduced.

In the case using Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

25 The distal tip of the guide wire is now positioned approximately 2 centimeters (cm) proximal from the distal tip of the nosecone. At this point the

application of vacuum as described in Step 2, Method 1 can be used to evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone.

Next, the cannula is advanced distally and guided out of the nosecone sideport via the internal exit ramps to pierce the sub-intimal tissue and gain access to the vessel true lumen. Once the vessel true lumen is accessed, the cannula remains in-place while the guide wire is advanced into the true lumen. The cannula may then be retracted into the catheter and resume its original position. The cannula may have various distal terminations, e.g. needle shaped, sharpened conical shaped, or circular serrations. The angle of the internal ramp can vary from approximately 30 to 80 degrees, but is not limited to these angles.

Lastly, the guide wire is held in position while the catheter is removed. As the catheter is retracted proximally over the guide wire, the floppy distal end of the guide wire may be able to pass through the nosecone side port. However, as the nosecone reaches the stiff mid and proximal sections of the guide wire, the guide wire falls through the slot connecting the side port with the end port. Therefore, as the catheter nosecone is retracted over the mid and proximal sections of the guide wire, it does so with the guide wire traveling through the nosecone distal port.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 to 0.060 inches; cannula element outside diameter is approximately 0.020 to 0.030 inches.

#### Embodiment 2 (Method 2 under Step 3) (Figure 8)

This embodiment includes: a nosecone or molded catheter termination attached to the distal end of the catheter lumen which houses the visualization element (Step 1, Methods 1 through 4); a lumen which houses either a guide wire

or the piercing element; lumens for vacuum ports; and optional forceps to physically secure the sub-intimal tissue.

The nosecone has a side port opening which allows for visualization of the vascular area (per Step 1, Methods 1 – 4) and identification of the orientation of the vessel true lumen. It also has an end port for tracking over a guide wire. Fluoroscopic features of the nosecone may also be used in conjunction with the primary visualization method to facilitate alignment to the true lumen. The nosecone also has a distal end port that allows the catheter to be tracked in a co-linear fashion over a conventional guide wire to the chosen vascular site.

The catheter shaft may be any of a number of catheter shafts known in the art. An example includes polymer lamination onto a stainless steel braided tube. Alternatively, the shafts may be fabricated as described in United States Patent Application Number 09/984,498, filed October 16, 2001. The nosecone may be fabricated/machined from stainless steel and EDM methods.

**Figure 31** is a flow diagram for piercing a pathway through sub-intimal tissue into a vessel true lumen, under an alternative embodiment. Procedurally, a guide wire is placed in the sub-intimal space of the target vasculature such that the guide wire distal end is located distal to the occluded area of the vessel. The pierce element is removed from the catheter, and the forceps (optional) are retracted to within the nosecone. The catheter is loaded onto the guide wire and tracked to the vascular site. The guide wire is removed. The pierce tool may be loaded into the catheter and the distal tip positioned just proximal to the nosecone sideport. The visualization element is activated and the sideport is directed towards the vessel true lumen.

Next, vacuum is applied to the vacuum lumens per Step 2, Method 2 thus beginning the process to invaginate the sub-intimal tissue into the nosecone

sideport. While confirming that the sideport remains directed towards the vessel true lumen, the skewers/forceps may be advanced thereby further securing the sub-intimal tissue per Step 2, Method 3.

The imaging element is then retracted, allowing more space for the sub-intimal tissue to further invaginate into the nosecone sideport. Note that under this configuration, visualization of the vessel true lumen is no longer possible since the imaging element is retracted proximal to the sideport.

Next, while maintaining vacuum, the pierce tool is advanced to pierce a pathway through the sub-intimal tissue, and through to the vessel true lumen. The pierce element is one of two fundamental types: including a lumen for a guide wire; and without a lumen. Each pierce element can have a variety of distal terminations, e.g., needle shaped, sharpened conical shaped or circular serrations around the outer diameter.

When using the pierce element without a lumen, once the pierce element has pierced a pathway across the sub-intimal tissue, it is retracted and removed from the catheter, and a guide wire introduced into the same catheter lumen and advanced through the sub-intimal tissue, and into the vessel true lumen.

When using the pierce element with an end lumen, once the pierce element has pierced a pathway across the sub-intimal tissue, its distal position is maintained within the vessel true lumen while a guide wire is advanced through the pierce element lumen. The guide wire exits the distal end of the pierce element and enters the vessel true lumen. The pierce element is then retracted fully into the nosecone.

For procedures using a pierce tool having a side lumen, after piercing, the guide wire is passed through the side port and into the vessel true lumen. The pierce tool is not retracted, and the catheter is removed over the guide wire.

Note that an alternative to the process described above is to not retract the visualization element prior to the piercing action. Thus the visualization element is positioned at the nosecone side port, and piercing is performed while simultaneously viewing the alignment to the vessel true lumen. Also note that in this alternative technique, the presence of the visualization element in the nosecone prevents the maximum invagination of tissue into the nosecone. Thus the first technique allows cutting through thicker sub-intimal tissue which separates the sub-intimal plane from the vessel true lumen.

Once a pathway is established into the vessel true lumen, the vacuum still applied at the nosecone aspirates blood from the vessel true lumen through the pathway in the sub-intimal tissue and into the catheter shaft, ultimately reaching the proximal end of the catheter. Further, the vacuum is lost. These two events may indicate that a pathway has been successfully established into the vessel true lumen. Once a pathway is established the guide wire is advanced into the true lumen, and the guide wire is held in position while the catheter is removed.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 to 0.060 inches; pierce element outside diameter is approximately 0.010 to 0.015 inches; and the imaging element outside diameter is as described in Step 1 above.

#### Method 3 under Step 3: Guidewire (Figures 4, and 7 through 15)

This method describes a catheter system which facilitates the use of either a conventional guide wire or specialized guide wire to establish a pathway across the sub-intimal tissue. The guide wire may be specially designed with a tip configuration that contains minute cutting members, or flutes, or the tip may be

processed to provide an abrasive surface. In either case, these tip surface features would cut or abrade a pathway through the sub-intimal tissue.

Micro-machining methods to fabricate the flutes or cutting features may include laser machining, electric discharge machining (EDM), or high-precision  
5 conventional machining. An abrasive tip surface may be fabricated using a micro-abrasive blaster using abrasive materials such as titanium oxide, or sodium bicarbonate. Very small abrasive features may also be laser machined on to the surface of the guide wire tip by the use of an Excimer laser in combination with a de-focusing mask. The de-focusing mask is a flat sheet fabricated from metal or  
10 other appropriate material, designed with a pattern of holes and/or slits or other shapes, which is placed between the laser light and the guide wire tip. This pattern is reduced in size and ablated onto the surface of the guide wire tip by the laser light that passes through the mask.

Note that these specialized guide wire tips are designed to cut or abraid a  
15 pathway through sub-intimal tissue when agitated/rotated and used in conjunction with the re-entry catheter, yet once introduced into the vessel true lumen, would not have the ability to exit into an extra-vascular space.

#### Embodiment 1 (Method 3 under Step 3) (Figure 4)

20 A first embodiment includes a distal nosecone and dual lumen catheter shaft. One lumen of the catheter shaft houses the imaging element, per Step 1, Methods 1-4, and the other lumen houses the guide wire.

Figure 32 is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under an embodiment.

25 Procedurally, the visualization element is removed from the catheter, and using

this lumen the catheter is tracked over a guide wire to the desired sub-intimal location.

Once the catheter is properly advanced in the sub-intimal plane, the guide wire is removed and the visualization element is advanced to the distal end of the nosecone. The visualization element is activated and the catheter is properly aligned to the vessel true lumen. Note that the pathway of the visualization element to the tissue may be through the shaft material itself. In the case of IVUS, this type of visualization may "see" through HDPE, and thus this is the preferred material for the dual lumen shaft for the visualization element lumen and the guide wire lumen. Alternatively, other visualization methods, e.g. Doppler, fiber optic, OCT may require a "window" through the visualization lumen to view the tissue.

Vacuum may be applied per Step 2, Method 1, evacuating the dissection plane and locking the sub-intimal tissue on the surface of the catheter. The guide wire is then pushed, and/or rotated to allow the guide wire distal tip to establish a pathway through the sub-intimal tissue and into the vessel true lumen. Lastly, while maintaining the guide wire position, the catheter is retracted and removed from the vasculature.

Typical dimensions of the catheter components are: outer shaft/nosecone outside diameter is approximately 0.050 to 0.070 inches; specialized guide wire outside diameter is approximately 0.010 to 0.018 inches; imaging element outside diameter is as described in Step 1 above.

#### Embodiment 2 (Method 3 under Step 3) (Figures 7 and 8)

This embodiment includes a distal nosecone and a multiple lumen catheter shaft that houses a visualization element per Step 1, Methods 1 through 4, a specialized guide wire, and optional separate vacuum ports.

Figure 33 is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a first alternative embodiment. Prior to the introduction of the catheter into the vasculature, the specialized guide wire is removed. Using this same lumen, the catheter is tracked  
5 over a guide wire to the desired sub-intimal location. Once the catheter is properly advanced in the sub-intimal plane, the guide wire is removed and replaced by the specialize guide wire.

The visualization element is activated and the catheter is properly aligned to the vessel true lumen. Vacuum may be applied per Step 2, Method 2, evacuating  
10 the dissection plane and invaginating the sub-intimal tissue into the catheter. Note that vacuum may be applied through the guide wire lumen, the visualization element lumen, or through optional vacuum ports as shown in Figure 7. At this point, the visualization element may be retracted proximally into the catheter shaft, adding more room for the sub-intimal tissue to be invaginated into the nosecone.  
15 This may be desired in the case that the sub-intimal tissue is thick, and requires a deeper purchase in order to create a pathway into the vessel true lumen. Note Figure 8 is a similar embodiment which shows optional forceps or skewers to hold the sub-intimal tissue.

Next, the guide wire is pushed, and/or rotated to allow the guide wire distal  
20 tip to establish a pathway through the sub-intimal tissue and into the vessel true lumen. While maintaining the guide wire position, the catheter is retracted and removed from the vasculature.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 to 0.060 inches;  
25 specialized guide wire outside diameter is approximately 0.010 to 0.018 inches; imaging element outside diameter is as described in Step 1 above.



**Embodiment 3 (Method 3 under Step 3) (Figure 9)**

This embodiment includes a nosecone or molded catheter termination attached to the distal end of the catheter, and an internal slidably disposed actuating cannula. The catheter shaft may be any of a number of catheter shafts known in the art. The nosecone includes a side exit port and a distal end port coupled via a slot which, and as will be described later, allows the guide wire to move from the side port into the distal end port when the catheter is retracted proximally over the guide wire. The cannula is guided out of the nosecone side port via internal exit ramps. The angle of the internal ramp is from approximately 30 degrees to 80 degrees, but not necessarily limited to these angles. The distal end port of the nosecone allows the catheter to be tracked in a co-linear fashion over a standard coronary guide wire. Representative dimensions are as follows: side port width is approximately 0.027 inches; slot width and distal port widths are approximately 0.016 inches; cannula outside diameter is approximately 0.025 inches; and cannula inside diameter is approximately 0.016 inches. Note that the internal ramp is the same width as the side port.

**Figure 34** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a second alternative embodiment. Procedurally, a guide wire is placed in the sub-intimal space of the target vasculature such that the guide wire distal end is located distal to the occluded area of the vessel. The cannula is retracted to a position proximal to the exit ramp so that the cannula exit port is co-linear with the inner diameter of the nosecone. This configuration allows the proximal end of the guide wire to be passed through the nosecone distal end port, the cannula and the catheter shaft, and thus the catheter may be tracked over the guide wire to the vascular site.

The catheter is then aligned to the vessel true lumen. This may be accomplished per Step 1 Methods 1-4, or Step 1, Method 5. In the case where Methods 1-4 are used, the guide wire (and optionally the cannula) is retracted from the catheter, and the visualization element is advance into the nosecone. The  
5 element is activated at the nosecone sideport and the side port is rotated to face the vessel true lumen. The visualization element is removed, and the cannula and guide wire are re-introduced. When Method 5 is used, the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

The distal tip of the guide wire is now positioned approximately 2  
10 centimeters proximal from the distal tip of the nosecone. At this point the application of vacuum as described in (Step 2, Method 1) may be used to evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone.

Next, the cannula is advanced distally and guided through the internal ramp  
15 until it is brought into secure purchase with the sub-intimal tissue. The guide wire is then advanced until the tip is coincident with the cannula distal tip, such that both are in contact with the sub-intimal tissue. Utilizing the combined effects of the vacuum and the slight extension of the cannula against the sub-intimal tissue, the guide wire is then rotated and pushed to initiate a pathway through the sub-intimal  
20 tissue.

After a re-entry pathway has been formed, the cannula is advanced into the true lumen, the RF system removed, and a conventional guide wire is placed into the true lumen. Alternatively, the cannula may remain retracted in the nosecone, and the guide wire fed directly through the pathway in the sub-intimal tissue and  
25 into the true lumen.

After the guide wire has been advanced into the vessel true lumen, the cannula may then be retracted into the catheter and resume its original position.

The position of the guide wire is maintained in the vessel true lumen, and the entire catheter system is retracted from the vasculature. As the catheter is retracted proximally over the guide wire, the floppy distal end of the guide wire may be able to pass through the nosecone side port, however as the nosecone reaches the stiff mid and proximal sections of the guide wire, the guide wire must now fall through the slot connecting the side port with the end port. Therefore, as the catheter nosecone is retracted over the mid- and proximal sections of the guide wire, it does so with the guidewire traveling through the nosecone distal port.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 to 0.060 inches; cannula element outside diameter is approximately 0.020 to 0.030 inches; and specialized guide wire outside diameter is approximately 0.010 to 0.018 inches.

#### Embodiment 4 (Method 3 under Step 3) (Figure 10)

This embodiment includes a catheter having a nosecone or molded distal termination, a single lumen catheter shaft, and a specialized guide wire to be used specifically with the catheter. This guide wire will be described in detail later. The catheter shaft may be any of a number of catheter shafts known in the art. The nosecone is similar to the nosecone of Figure 9 in that it includes a side exit port and a distal end port coupled via a slot. However, the dimensions of these features are peculiar to this design.

As previously stated, the nosecone of Figure 10 is similar to that of Figure 9 with the exception that the width of the side port and the outside diameter of the end port are slightly larger than the width of the connecting slot. As example

dimensions, the width of the side port and the outside diameter of the end port may be approximately 0.016 inches, and the slot may be approximately 0.012 inches. The significance of these dimensions will become evident once the guide wire is dimensionally described. The internal ramp has a width approximately  
5 equal to the size of the side port. The angle of the internal ramp may also vary from approximately 30 degrees to 80 degrees, but is not necessarily limited to these angles. The distal end port of the nosecone allows the catheter to be tracked in a co-linear fashion over the specialized guide wire.

The guide wire may be fabricated using standard methods and materials  
10 known by those skilled in the art. The outside diameter over the distal most 5-8 centimeters may be approximately 0.014 inches, followed proximally by a 1-2 centimeter portion with an outside diameter of approximately 0.010 inches, followed proximally by the remainder of the guide wire with an outside diameter of approximately 0.014 inches.

15 **Figure 35** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a third alternative embodiment. Procedurally, the catheter is tracked over a standard guide wire, or over the specialized guide wire, either of which has been placed in the desired sub-intimal space of the target vasculature. If a standard guide wire is initially  
20 used, it may then be removed and the specialized guide wire advanced into the catheter.

Next, the catheter must be aligned to the vessel true lumen. This may be accomplished per Step 1 Methods 1-4, or Step 1, Method 5. In the case where Methods 1-4 are used, the guide wire is retracted from the catheter, and the  
25 visualization element is advance into the nosecone. The element is activated at the nosecone sideport and the side port is rotated to face the vessel true lumen.

The visualization element is removed, and the specialized guide wire is re-introduced. In the case that Method 5 is used, the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

Next, the distal tip of the specialized guide wire is positioned approximately  
5 2 centimeters proximal from the distal tip of the nosecone. At this point the application of vacuum as described in (Step 2, Method 1) may be used to evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone.

The specialized guide wire is now advanced distally. Note that because the  
10 nosecone exit port is co-incident with the catheter lumen, the natural tendency of the specialized guide wire may be for it to simply pass out the distal port of the nosecone. The intent, however is to advance the specialized guide wire out of the nosecone side port. This issue is easily resolved. Prior to the introduction of any guide wire into the vasculature, the physician routinely places a small curve on the  
15 end of the guide wire to facilitate negotiating the tortuosity in the vasculature. Therefore, the specialized wire need only be rotated until the curved distal portion of the wire engages the internal rails. It may then be advanced onto the internal ramp of the side port.

The distal 5-8 centimeters of the specialized guide, at a width of  
20 approximately 0.014 inches, will advance on the two rails of the internal ramp, since the rails are separated by the 0.012-inch width of the slot, and the side port is 0.016 inches wide. The specialized guide wire is brought into contact with the sub-intimal tissue. Utilizing the effect of the vacuum, the specialized guide wire can be rotated and pushed to initiate a pathway through the sub-intimal tissue.

25 After the specialized guide wire has crossed to the vessel true lumen, the wire is advanced further distally. In this process, the 5-8 centimeters of distal

length guide wire (0.014 inches outside diameter) continues to advance over the internal rails, until the 1 centimeter section of guide wire (0.010 inches outside diameter) reaches the distal edge of the ramp. At this point the section of guide wire having a 0.010 inch outside diameter falls through the 0.012-inch wide slot  
5 and into the 0.016-inch wide distal end port.

At this point, the distal 6-9 centimeters of the guide wire is across the sub-intimal plane and into the vessel true lumen. While maintaining the wire position in the vasculature, the catheter may now be retracted proximally along the wire, because the nosecone end port is 0.016 inches in diameter, and the proximal  
10 portion of the wire is 0.014 inches in diameter.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; specialized guide wire outside diameter is as described above.

15        Embodiment 5 (Method 3 under Step 3) (Figure 11)

This embodiment includes a catheter having a simple nosecone or molded distal termination, and a single lumen catheter shaft. The catheter shaft may be any of a number of catheter shafts known in the art. The nosecone includes an internal ramp connecting the catheter lumen with a single side exit port. This  
20 embodiment has no distal port to track over a guide wire. This catheter can be used in conjunction with a conventional guide wire, or a specialized guide wire as described at the end of this section.

Figure 36 is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a fourth alternative  
25 embodiment. Procedurally, the catheter is tracked over a standard guide wire which has been placed in the desired sub-intimal space of the target vasculature.

Note that since the guide wire emerges laterally from the nosecone, the very tip of the catheter will track eccentrically over the guide wire.

Next, the catheter is aligned to the vessel true lumen. This may be accomplished per Step 1 Methods 1-4, or Step 1, Method 5. Figure 11 shows both configurations. In the case of Methods 1-4, the guide wire is retracted from the catheter, and the visualization element is advanced into the nosecone. The element is activated at the nosecone side port and the side port is rotated to face the vessel true lumen. The visualization element is removed, and the guide wire is re-introduced. In the case of Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

Next, the distal tip of the guide wire is positioned approximately 2 centimeters proximal from the distal tip of the nosecone. At this point the application of vacuum as described in (Step 2, Method 1) may be used to evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone. Vacuum is translated to the nosecone via the single shaft lumen.

The guide wire is then advanced distally into the nosecone. The guide wire is brought into contact with the sub-intimal tissue at the nosecone side port at an angle determined by the exit ramp of the nosecone. The wire is then rotated and pushed in order to initiate and propagate a pathway through the sub-intimal tissue. The angle of the internal ramp may vary from approximately 30 degrees to 80 degrees, but is not necessarily limited to these angles.

After the guide wire has successfully been advanced into the vessel true lumen, the guide wire position is maintained and the catheter may be retracted proximally over the guide wire and removed from the vasculature.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; specialized guide wire outside diameter is approximately 0.010 to 0.018 inches.

5            Embodiment 6 (Method 3 under Step 3) (Figure 12)

     This embodiment includes a catheter having a nosecone or molded distal termination, a single lumen catheter shaft, and an internal slidably disposed tube the distal end of which translates within the catheter nosecone. The distal end of the push tube or member is slidably disposed within the nosecone. Upon  
10    actuation of the push member in a distal direction, a percentage of the proximal section of the nosecone side port is covered. The percentage of coverage is controlled by a distal stop within the nosecone that limits the distal translation of the internal sliding member or tube. In the case the internal sliding member is a tube, the internal sliding member becomes the guide wire lumen.

15            As the internal sliding member is advanced into this distal position, it reduces the effective length of the nosecone side port. This forces the guide wire to exit the side port at a more acute angle that is more normal to the axis of the catheter, and more normal to the sub-intimal tissue plane. The exit angle of the guide wire is governed by the proximal contact point against the internal sliding  
20    member, and the distal contact point against the exit ramp of the nosecone. The greater acute angle allows the guide wire to produce more force normal to the sub-intimal tissue surface, and improves the ability of the wire to establish a pathway across the sub-intimal tissue.

**Figure 37** is a flow diagram for using a guide wire to establish a pathway  
25    through sub-intimal tissue into a vessel true lumen, under a fifth alternative embodiment. Procedurally, the sliding tube may remain in the catheter at all times.



In preparation of the catheter and during delivery of the catheter to the vascular site, the tube is retracted just proximal to the nosecone. Prior to the advancement of the guide wire out of the nosecone sideport, the sliding tube is advanced to its distal most position, thus reducing the effective length of the nosecone sideport.

5        Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; internal slide tube outside diameter is approximately 0.020 to 0.030 inches; and specialized guide wire outside diameter is approximately 0.010 to 0.018 inches.

10        Embodiment 7 (Method 3 under Step 3) (Figure 13)

Figure 13 is a diagram of a single lumen catheter shaft, terminated in a "J" tip, used in conjunction with a conventional or specialized guide wire. The "J" tip configuration of the catheter is designed to be torqued into position within the sub-intimal plane and directed towards the sub-intimal tissue. The guide wire is  
15        directed at an angle normal to the sub-intimal tissue, which improves the ability of the wire to establish a pathway across the sub-intimal tissue.

The "J" termination of the catheter can be fluoroscopically visible to facilitate the positioning process in the sub-intimal plane. This type of termination may be easily fabricated or molded from fluoroscopically visible materials such as platinum  
20        coils or gold coated stainless steel coils laminated with a variety of polymers, e.g. nylons, HDPE, or Pebax. The "J" tip is designed to straighten in order to track over a guide wire to the vascular site, yet re-form its shape when positioned in the vasculature, and the guide wire is retracted. A visualization window may be incorporated just proximal to the "J" tip to be used in conjunction with an on-board  
25        visualization technique, per Step 1, Methods 1-4. For the use of IVUS, for example, this window may be fabricated from HDPE.

**Figure 38** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a sixth alternative embodiment. Procedurally, a guide wire is positioned in the sub-intimal space. The distal end of the catheter is loaded onto the guide wire. In this process the "J" tip is straightened as it tracks over the wire and to the sub-intimal site. Once the terminal end of the catheter has reached the desired location, the guide wire is retracted, allowing the "J" tip to re-form.

The catheter is now aligned to the vessel true lumen. This may be accomplished per Step 1 Methods 1-4, or Step 1, Method 5. Figure 13 shows both configurations. In the case of Methods 1-4, the guide wire is retracted from the catheter, and the visualization element is advanced to the visualization window. The visualization element is activated at the window and the "J" tip is rotated to face the vessel true lumen. The visualization element is removed, and the guide wire is re-introduced. In the case of Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

Next, vacuum may be applied through the catheter lumen per Step 2, method 1 to evacuate the sub-intimal plane. The guide wire is then advanced until it is brought into contact with the sub-intimal tissue. Procedurally, the wire is then pushed or rotated as required in order to initiate and propagate a pathway through the sub-intimal tissue and into the vessel true lumen. Once the guide wire has established a pathway through the sub-intimal tissue, the wire position is maintained, and the catheter is retracted, leaving the distal portion of the wire positioned in the vessel true lumen.

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; specialized guide wire outside diameter is approximately 0.010 to 0.018 inches.

**Embodiment 8 (Method 3 under Step 3) (Figure 14)**

This embodiment includes a catheter having a nosecone or molded distal termination, an internal push-ramp which is actuated by an internal push tube or member, and a single lumen catheter shaft. The push ramp may be constructed of  
5 a flexible metal such as Nitinol or spring steel, or a polymer such as nylon or PEEK, any of which fabricated with or without appropriate detents to allow for bending, as required. The distal end of the push ramp is connected or hinged in some fashion about the internal distal termination of the catheter shaft, opposite the side port. When the push tube or member is fully retracted proximally, the  
10 push ramp assumes a linear configuration, lying essentially flat against the inside wall of the catheter, opposite the nosecone side port. When the push tube or member is advanced distally, the push ramp forms an incline that leads from the proximal end of the ramp, opposite the nosecone side port, to the distal end of the side port. This ramp will re-direct the guide wire out of the nosecone side port as it  
15 is advanced distally in the catheter.

This catheter system may be used with any of the visualization techniques of Step 1.

**Figure 39** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a seventh alternative  
20 embodiment. Procedurally, a guide wire is positioned in the sub-intimal space. The pull tube or member is retracted proximally, and the distal end of the catheter is loaded onto the guide wire. Once the terminal end of the catheter has reached the desired location, the guide wire is retracted just proximal to the proximal portion of the ramp.

25 The catheter is aligned to the vessel true lumen. This may be accomplished per Step 1 Methods 1-4, or Step 1, Method 5. Figure 13 shows both

configurations. In the case of Methods 1-4, the guide wire is retracted from the catheter, and the visualization element is advanced to the visualization window.

The visualization element is activated at the nosecone sideport and the sideport is rotated to face the vessel true lumen. The visualization element is removed, and

5 the guide wire is re-introduced. In the case of Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

— Next, vacuum may be applied through the catheter lumen per Step 2, method 1 to evacuate the sub-intimal plane. The push tube or member is then advanced distally to urge the push ramp into its hinged configuration. The guide  
10 wire is then advanced distally, following the ramp to the nosecone cone port until it is brought into contact with the sub-intimal tissue. Procedurally, the wire is then pushed or rotated as required in order to initiate and propagate a pathway through the sub-intimal tissue and into the vessel true lumen. Once the guide wire has established a pathway through the sub-intimal tissue, the wire position is  
15 maintained and the catheter is retracted, leaving the distal portion of the wire positioned in the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; internal push tube outside diameter is approximately 0.020 to 0.030 inches; and  
20 specialized guide wire outside diameter is approximately 0.010 to 0.018 inches.

#### Embodiment 9 (Method 3 under Step 3) (Figure 15)

This embodiment includes a dual lumen catheter shaft terminating in a single lumen "J" tip, and is used in conjunction with a conventional or specialized  
25 guide wire. This dual lumen catheter has the same "J" type single lumen distal termination as described in Method 3, Embodiment 7, with the exception that it

transitions proximally to a dual lumen catheter shaft. Fabrication and materials for the "J" tip are as stated in Method 3, Embodiment 7. The dual lumen catheter shaft may be fabricated using standard materials and methods known to those skilled in the art. Only one of the slidably disposed elements contained within  
5 either lumen can be advanced individually into the "J" tip single lumen, as required by the procedure. The "J" tip may accommodate only one element at any time.

For example, one lumen may contain the guide wire while the other lumen may be used to deliver various elements to the vascular site, e.g. visualization elements per Step 1, Methods 1-4, or other types of re-entry elements such as a  
10 wire with a stiff distal tip which could be used only to pierce a hole through the sub-intimal tissue, but would be too stiff to be advanced into the vessel true lumen.

**Figure 40** is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under an eighth alternative embodiment. In a first scenario, the first lumen of the dual lumen contains a  
15 standard guide wire and the second lumen contains a re-entry wire or re-entry element to pierce or otherwise establish a pathway into the vessel true lumen. Procedurally, one lumen is loaded with a re-entry wire or element and advanced just proximal to the entrance to the single distal lumen.

The distal end of the catheter is then loaded onto a standard guide wire  
20 which has been advanced into a sub-intimal plane. The catheter is advanced to the desired vascular location, and the guide wire withdrawn just proximal to the entrance to the single lumen. The catheter is now aligned to the vessel true lumen. This embodiment would make use of Step 1, Method 5 to align the catheter.

25 The sub-intimal plane is now evacuated per Step 2, Method 2. Next, the guide wire or re-entry element is advanced into the "J" tip and manipulated to

establish a pathway through the sub-intimal plane and into the vessel true lumen.

The re-wire or re-entry element is retracted out of the single lumen. The standard guide wire may then be advanced into the distal single lumen, and out of the "J" tip, through the pathway produced in the sub-intimal tissue, and into the vessel true lumen. Lastly, the guide wire is held in place while the catheter is retracted proximally and removed from the vasculature.

Figure 41 is a flow diagram for using a guide wire to establish a pathway through sub-intimal tissue into a vessel true lumen, under a ninth alternative embodiment. In a second scenario, the first lumen of the dual lumen contains a visualization element, per Step 1, Method 1-4, and the second lumen contains a re-entry wire to establish a pathway into the vessel true lumen.

Procedurally, one lumen is loaded with the visualization element and advanced proximal to the entrance to the distal single lumen. The distal end of the catheter is then loaded onto a standard guide wire that has been advanced into a sub-intimal plane. The catheter is advanced to the desired vascular location, the guide wire is removed from the catheter and replaced with the re-entry wire. Next the visualization element is advanced into the distal single lumen within the area of the window for viewing, and the "J" tip is aligned with the vessel true lumen. The visualization element is then withdrawn to within the dual lumen portion of the shaft.

The sub-intimal plane is now evacuated per Step 2, Method 2. This is best accomplished through the lumen which houses the visualization element because the visualization element need not be translated proximally/distally through the proximal hemostasis seal for the remainder of the procedure. Next, the re-entry wire is advanced into the "J" tip and manipulated to establish a pathway through the sub-intimal tissue and into the vessel true lumen. The re-entry wire is retracted

and replaced with the guide wire, and the guide wire is held in place while the catheter is retracted proximally and removed from the vasculature.

In a third scenario, all steps are similar to the second scenario, with the exception that a stiff re-entry wire or element is used to establish the pathway  
5 through the sub-intimal tissue, and is then replaced by a standard guide wire and advanced into the vessel true lumen before removal of the catheter.

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; dual lumen shaft is approximately 0.030 to 0.050 inches (each lumen); and specialized guide wire  
10 outside diameter is approximately 0.010 to 0.018 inches.

**Method 4 under Step 3: Radio Frequency (RF) Energy (Figures 5, 6, 7, 9, 11 through 15, and 18)**

This method describes the application of radio frequency (RF) energy to  
15 ablate a select section of sub-intimal tissue to create a path into the true lumen. The RF energy can be delivered as continuous with an unspecified duration. Alternatively, the RF energy may be gated such that it is defined by a pre-determined duration, e.g., 5-50 milliseconds.

Two separate modes of application of RF energy may be employed,  
20 unipolar or bi-polar. Both methods may be employed in all embodiments described herein.

In a unipolar configuration, a single active electrode (or a group of common electrodes) is located at a distal position of the re-entry catheter. The active electrode(s) may be contained on the outside surface of the catheter, or within the  
25 distal nosecone or distal termination of the catheter, such that a communication path exists to the outside of the catheter. The placement of the active electrode(s)

is such that when the catheter is aligned to the vessel true lumen per Step 1, the electrode(s) faces the sub-intimal tissue to be ablated. A second grounding electrode is external to the patient and placed against the patient's buttocks. The surface area of the grounding electrode is very large with respect to the catheter active electrode.

As RF energy is applied between the catheter active electrode(s) and the grounding electrode, a closed circuit is established from one electrode to the other, through the patient's body tissue. The energy density (RF power) and duration is adjusted to ablate only the sub-intimal tissue that separates the dissection plane from the vessel true lumen. As the RF energy travels between the catheter electrode(s) and the grounding electrode, the energy density decreases significantly such that immediately peripheral to the volume of sub-intimal tissue to be ablated, the energy density is insufficient to affect the surrounding vessel structure or other body tissues.

In a bipolar configuration both the active electrode and the grounding electrode are contained at the distal portion of the catheter. One or both can be mounted on the outside surface of the catheter, or within the catheter itself. Both electrodes are intended to be of similar size, although exact size and configuration may vary to meet the overall design requirements of the re-entry catheter system. The placement of the electrodes is such that when the catheter is aligned to the vessel true lumen per Step 1, the sub-intimal tissue to be ablated completes the circuit between electrodes.

As RF energy is applied between the catheter electrodes, a closed circuit is established from one electrode to the other, co-linear with the sub-intimal tissue plane that separates the dissection plane from the vessel true lumen. The energy density (RF power) and duration is adjusted to ablate the sub-intimal tissue that



separates the dissection plane from the vessel true lumen. Unlike the unipolar technique described above, the energy density between the two bipolar electrodes is relatively constant, and only the sub-intimal tissue along the path between the two electrodes is ablated. It is surmised that the bipolar mode may have more accurate control over tissue ablation than the unipolar configuration.

Embodiment 1 (Method 4 under Step 3) (Figures 5 through 7A)

This embodiment includes a catheter having a nosecone or molded distal termination, an internal visualization element, and one or more RF electrodes. In this embodiment the single lumen catheter shaft is terminated by a formed nosecone having a sideport for imaging the sub-intimal tissue and locating the vessel true lumen, and an endport for tracking the catheter over a guide wire. The nosecone is also shown with one or more RF electrodes, suitable to embody either the unipolar or bipolar configuration.

**Figure 42** is a flow diagram for creating a path into a vessel true lumen using radio frequency (RF) energy, under an embodiment. Procedurally, the visualization element is removed from the catheter (Fig. 5 and 6 only), and the catheter is loaded onto a guide wire that has been advanced into a sub-intimal plane. The catheter is tracked to the vascular site, and the guide wire is removed. The visualization element is loaded into the catheter, and advanced into the nosecone sideport. The visualization element is activated and the sideport is directed towards the vessel true lumen.

Next, vacuum is applied within the catheter lumen per Step1, Method 1 to secure the sub-intimal tissue to the nosecone, and bring the sub-intimal tissue into contact with the RF electrodes. Next, RF energy may be applied to the electrodes

as required to ablate the sub-intimal tissue and form a pathway into the vessel true lumen.

If the sub-intimal tissue is thick, it may be advantageous to invaginate more sub-intimal tissue into the nosecone by retracting the visualization element and continuously applying vacuum per Step 2, Method 2. This enhances the likelihood of ablating the volume of sub-intimal tissue necessary to establish a pathway into the true lumen.

As addressed in the beginning of this Method, RF energy may be applied in the unipolar configuration to one or more of the electrodes mounted on the catheter, treating the electrodes as a common pole with respect to the grounding plate. Alternatively, in the bipolar configuration, the RF signal may be applied to one of the catheter mounted electrodes, using the other as the ground return.

The visualization element is removed (Fig. 5 and 6 only), and a standard guide wire is introduced either through the nosecone end port or side port into the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.040 to 0.060 inches; imaging element outside diameter is as described in Step 1 above.

#### **Embodiment 2 (Method 4 under Step 3) (Figure 9)**

This embodiment includes a nosecone or molded catheter termination attached to the distal end of the catheter, an internal slidably disposed actuating cannula, and an element slidably disposed in the cannula which contains distal electrodes.

**Figure 43** is a flow diagram for creating a path into a vessel true lumen using radio frequency (RF) energy, under a first alternative embodiment. This

embodiment is essentially similar to that of Method 3, Embodiment 3 with the exception that in place of the guide wire, an element is delivered inside the cannula which contains one or more electrodes to ablate the sub-intimal tissue.

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; cannula element outside diameter is approximately 0.020 to 0.030 inches; and RF element outside diameter is approximately 0.015 to 0.020 inches.

#### Embodiment 3 (Method 4 under Step 3) (Figure 11)

This embodiment includes a catheter having a simple nosecone or molded distal termination, a single lumen catheter shaft and an element slidably disposed in the catheter which contains distal electrodes.

Figure 44 is a flow diagram for creating a path into a vessel true lumen using radio frequency (RF) energy, under a second alternative embodiment.

Procedurally, this embodiment is essentially similar to that of Method 3, Embodiment 5 with the exception that in place of the guide wire, an element is delivered inside the catheter which contains one or more distal electrodes to ablate the sub-intimal tissue.

After a pathway has been established through the sub-intimal tissue, the slidably disposed element is removed from the catheter, and a standard guide wire is advanced into the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; and RF element outside diameter is approximately 0.015 to 0.020 inches.

**Embodiment 4 (Method 4 under Step 3) (Figure 12)**

This embodiment includes a catheter having a simple nosecone or molded distal termination, a single lumen catheter shaft, an internal slidably disposed tube, and an element slidably disposed in the tube which contains distal electrodes.

5       **Figure 45** is a flow diagram for creating a path into a vessel true lumen using radio frequency (RF) energy, under a third alternative embodiment. Procedurally, this embodiment is essentially similar to that of Method 3, Embodiment 6, with the exception that in place of the guide wire, a slidably disposed element is used which contains one or more electrodes to ablate the  
10   sub-intimal tissue.

After a pathway has been established through the sub-intimal tissue, the slidably disposed element is removed from the catheter, and a standard guide wire is advanced into the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer  
15   shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; internal slide tube outside diameter is approximately 0.020 to 0.030 inches; and RF element outside diameter is approximately 0.015 to 0.020 inches.

---

**Embodiment 5 (Method 4 under Step 3) (Figure 13)**

20       This embodiment includes a single lumen catheter shaft, terminated in a "J" tip, used in conjunction with an element slidably disposed in the catheter which contains distal electrodes. Procedurally, this embodiment is used in a similar fashion to Method 3, Embodiment 7. In this embodiment, and with further reference to **Figure 38**, the slidably disposed element contains one or more  
25   electrodes to ablate the sub-intimal tissue, and takes the place of the guide wire.

After a pathway has been established through the sub-intimal tissue, the slidably disposed element is removed from the catheter, and a standard guide wire is advanced into the vessel true lumen

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; and RF element outside diameter is approximately 0.015 to 0.020 inches.

---

**Embodiment 6 (Method 4 under Step 3) (Figure 14)**

This embodiment includes a catheter having a nosecone or molded distal termination, an internal push-ramp which is actuated by an internal push tube or member, an element slidably disposed in the catheter which contains distal electrodes and a single lumen catheter shaft.

Procedurally, this embodiment is used in a similar fashion to Method 3, Embodiment 8. In this embodiment, and with further reference to **Figure 39**, the slidably disposed element contains one or more electrodes to ablate the sub-intimal tissue, and takes the place of the guide wire.

After a pathway has been established through the sub-intimal tissue, the slidably disposed element is removed from the catheter, and a standard guide wire is advanced into the vessel true lumen

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; internal push tube outside diameter is approximately 0.020 to 0.030 inches; and RF element outside diameter is approximately 0.015 to 0.020 inches.

**Embodiment 7 (Method 4 under Step 3) (Figure 15)**

Describes a dual lumen catheter shaft, terminated in a single lumen "J" tip, used in conjunction with an element slidably disposed in the catheter which contains distal electrodes. Procedurally, this embodiment is used in a similar fashion to Method 3, Embodiment 9. In this embodiment, and with further reference to **Figures 40 and 41**, the slidably disposed element contains one or more electrodes to ablate the sub-intimal tissue, and takes the place of the guide wire.

After a pathway has been established through the sub-intimal tissue, the slidably disposed element is removed from the catheter, and a standard guide wire is advanced into the vessel true lumen

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; dual lumen shaft outside diameter is approximately 0.030 to 0.050 inches (each lumen); and RF element outside diameter is approximately 0.015 to 0.020 inches.

**Method 5 under Step 3: Laser Energy (Figures 4, 7 through 15; and 17)**

This method describes the application of laser energy to ablate the sub-intimal tissue that separates the dissection plane from the vessel true lumen. In the following embodiments, optical fiber(s) is/are contained within the catheter, such that the terminal end of the optical fiber(s) is/are positioned to deliver laser energy to ablate the sub-intimal tissue.

In the following embodiments, the optical fiber element may comprise: a single optical fiber contained within a protective sheath such as polyimide or high density polyethylene; a bundle of optical fibers collectively protected by a sheath such as polyimide, or high density polyethylene; or an optical fiber bundle

arranged in an annular fashion, such that the interior of the bundle is an open lumen sized to accommodate a standard guide wire. In this embodiment, the optical fibers may be arranged in the annular fashion within a polymer extrusion. The extrusion design thus encases the fibers in the annular arrangement, provides  
5 a lumen within the extrusion for the passage of a guide wire, and provides an outer smooth surface, such that the fiber optic bundle may be translated within a catheter lumen. These optical fiber embodiments are shown in Figure 16.

Embodiment 1 (Method 5 under Step 3) (Figure 4)

10 This embodiment includes a distal nosecone or molded shaft termination and dual lumen catheter shaft. A first lumen of the catheter shaft houses the imaging element, per Step 1, Method 1-4, and a second lumen houses the optical fiber system or a guide wire.

This embodiment allows simultaneous visualization of the true lumen and  
15 sub-intimal tissue, while delivering laser energy via the optical fiber system to ablate the sub-intimal tissue. In operation, small volumes of sub-intimal tissue are sequentially ablated with each pulsed delivery of optical energy, until a pathway has been established. The exit angle from the nosecone of the lumen which houses the optical fiber system is approximately in the range of 20 degrees to 90  
20 degrees, but is not necessarily limited to these angles.

The termination of the optical fiber may be normal to the axis of the optical fiber, such that the exit angle is zero degrees with respect to the optical fiber. Alternatively, the termination of the optical fiber may be angled with respect to the axis to provide total internal reflection of the light, such that the light emerges at an  
25 angle to the axis of the optical fiber. The angle termination provides a greater overall exit angle of the laser light, when combined with the angle by which the

optical fiber exits the corresponding lumen, and provides overall operational flexibility.

Figure 46 is a flow diagram for creating a path into a vessel true lumen using laser energy, under an embodiment. Procedurally, prior to introducing the catheter into the vasculature, the optical fiber system is loaded into a lumen, and the visualization element is removed from the catheter. Using the visualization lumen, the catheter is loaded onto a guide wire that has been advanced into a sub-intimal plane, and the catheter is advanced to the desired vascular site. In this configuration, the catheter tracks on the wire in a co-linear fashion.

10 The guide wire is next removed from the catheter and the visualization element is advanced in the same lumen until it is properly positioned at the distal end of the catheter. The optical fiber system is advanced until the distal termination is coincident with the lateral exit port at the nosecone.

Next, the visualization element is activated and the catheter is properly aligned to the vessel true lumen. Note that the pathway from the visualization element to the tissue may be through the shaft material itself. In the case of IVUS, this type of visualization may "see" through HDPE, and thus is the preferred material for the dual lumen shaft, the visualization element lumen, and the guide wire lumen. Alternatively, other visualization methods, e.g., Doppler, fiber optic, and OCT may require a "window" from the visualization lumen to view the tissue.

20 Vacuum is now applied per Step 2, Method 1, evacuating the dissection plane and locking the sub-intimal tissue onto the surface of the catheter. Vacuum may be applied through the visualization element lumen, the optical fiber lumen, or both lumens.

25 Laser energy is delivered to the optical fiber system as required to ablate the sub-intimal tissue which separates the dissection plane from the vessel true



lumen. If the fiber optic system that does not incorporate the guide wire is used, the fiber optic system is removed and a standard guide wire is introduced for advancement through the pathway in the sub-intimal tissue and into the true lumen. In the case of the fiber optic system that contains a guide wire lumen, a  
5 guide wire is simply advanced through the lumen, through the pathway in the sub-intimal tissue, and into the true lumen.

Typical dimensions of the catheter-components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 to 0.070 inches; fiber optic system outside diameter is approximately 0.010 to 0.030 inches; and imaging  
10 element outside diameter is as described in Step 1 above.

**Embodiment 2 (Method 5 under Step 3) (Figures 7 and 8)**

The embodiment of **Figure 7** includes a distal nosecone and a multi-lumen catheter shaft which house a visualization element per Methods 1-4 under Step 1,  
15 an optical fiber system, and optional separate vacuum ports. The embodiment of **Figure 8** includes optional sub-intimal tissue forceps.

**Figure 47** is a flow diagram for creating a path into a vessel true lumen using laser energy, under a first alternative embodiment. Prior to the introduction of the catheter into the vasculature, the optical fiber system is removed from a  
20 corresponding lumen. The catheter is tracked over a guide wire within this lumen to the desired sub-intimal location. Once the catheter is properly advanced in the sub-intimal plane, the guide wire is removed and the optical fiber system is replaced.

The visualization element is now activated and the catheter is properly  
25 aligned to the vessel true lumen. Next, vacuum may be applied per Step 2, Method 2, thereby evacuating the dissection plane and invaginating the sub-intimal

tissue into the catheter nosecone or distal termination. Note that vacuum may be applied through the fiber optic lumen, the visualization element lumen, or through optional vacuum ports like those shown in Figures 7 and 8.

While maintaining alignment to the vessel true lumen, the visualization  
5 element may be retracted proximally into the catheter shaft, making more room for the sub-intimal tissue to be invaginated into the nosecone. This may be desired in the case that the sub-intimal tissue is thick, and requires a deeper purchase in order to create a pathway into the vessel true lumen. Figure 8 is a similar embodiment showing optional forceps or skewers holding the sub-intimal tissue.

10 Next, laser energy is delivered to the optical fiber system as required to ablate the sub-intimal tissue that separates the dissection plane from the vessel true lumen. If a fiber optic system that does not incorporate the guide wire is used, the fiber optic system is removed and a standard guide wire is introduced for advancement through the pathway in the sub-intimal tissue and into the true  
15 lumen. In the case of the fiber optic system that contains a guide wire lumen, a guide wire is simply advanced through the lumen, through the pathway in the sub-intimal tissue, and into the true lumen.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 to 0.060 inches; fiber  
20 optic system outside diameter is approximately 0.010 to 0.030 inches; and imaging element outside diameter is as described in Step 1 above.

### Embodiment 3 (Method 5 under Step 3) (Figure 9)

This embodiment includes a nosecone or molded catheter termination  
25 attached to the distal end of the catheter, an internal slidably disposed actuating cannula, and an optical fiber system.

The catheter shaft may be any of a number of catheter shafts known in the art. The nosecone includes a side exit port and a distal end port coupled via a slot which, and as is described below, allows the guide wire to translate or move from the side port into the distal end port when the catheter is retracted proximally over the guide wire. The cannula is guided out of the nosecone sideport via an internal exit ramp. The cannula tip of an embodiment is bluntly terminated, but is not so limited. The angle of the internal ramp may vary from approximately 30 degrees to 80 degrees, but is not necessarily limited to these angles. The distal end port of the nosecone allows the catheter to be tracked in a co-linear fashion over a standard coronary guide wire. Representative component dimensions of an embodiment are as follows: side port width is approximately 0.027 inches; slot width and distal port width are both approximately 0.016 inches; and cannula outside diameter is approximately 0.025 inches and inside diameter is approximately 0.016 inches. Note that the internal ramp is the same width as the side port.

**Figure 48** is a flow diagram for creating a path into a vessel true lumen using laser energy, under a second alternative embodiment. Procedurally, a guide wire is placed in the sub-intimal space of the target vasculature such that the guide wire distal end is located distal to the occluded area of the vessel. The cannula is retracted to a position proximal to the exit ramp so that the cannula exit port is co-linear with the inner diameter of the nosecone. This configuration allows the proximal end of the guide wire to be passed through the nosecone distal end port, the cannula and the catheter shaft. Thus the catheter may be tracked over the guide wire to the vascular site.

Next, the catheter is aligned to the vessel true lumen. This may be accomplished per under Step 1, Methods 1-4 or Step 1, Method 5. In the scenario

where any of Methods 1-4 are used, the guide wire (and optionally the cannula) is retracted from the catheter, and the visualization element is advanced into the nosecone. The element is activated at the nosecone side port and the side port is rotated to face the vessel true lumen. The visualization element is then removed.

- 5 When using Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization. The guide wire is then removed.

— The cannula and fiber optic system are re-introduced. The distal tip of the optical fiber system is positioned approximately 2 centimeters proximal from the distal tip of the nosecone. At this point, the application of vacuum can be used to  
10 evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone.

Next, the cannula is advanced distally and guided through the internal ramp until it is brought into secure purchase with the sub-intimal tissue. The optical fiber system is then advanced until the tip is coincident with the cannula distal tip such  
15 that both are in contact with the sub-intimal tissue.

Laser energy is now delivered to the optical fiber system as required to ablate the sub-intimal tissue that separates the dissection plane from the vessel true lumen. If the fiber optic system that does not incorporate the guide wire is used, the fiber optic system is removed and a standard guide wire is introduced for  
20 advancement through the pathway in the sub-intimal tissue and into the true lumen. In the case of the fiber optic system that contains a guide wire lumen, a guide wire is advanced through the lumen, through the pathway in the sub-intimal tissue, and into the true lumen. The cannula is subsequently retracted to its original position, so that its distal end is positioned just proximal to the nosecone  
25 internal ramp.

The position of the guide wire is maintained in the vessel true lumen, and the entire catheter system is retracted from the vasculature. As the catheter is retracted proximally over the guide wire, the floppy distal end of the guide wire may be able to pass through the nosecone side port. However as the nosecone reaches the stiff mid and proximal sections of the guide wire, the guide wire falls through the slot connecting the side port with the end port. Therefore, as the catheter nosecone is retracted over the mid- and proximal sections of the guide wire, it does so with the guide wire traveling through the nosecone distal port.

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; cannula element outside diameter is approximately 0.020 to 0.030 inches; and fiber optic system outside diameter is approximately 0.010 to 0.030 inches.

#### Embodiment 4 (Method 5 under Step 3) (Figure 11)

This embodiment includes a catheter having a simple nosecone or molded distal termination, a single lumen catheter shaft, and a fiber optic system slidably disposed in the catheter lumen. The catheter shaft can be any of a number of catheter shafts known in the art. The nosecone includes an internal ramp connecting the catheter lumen with a single side exit port. This design has no distal port to track over a guide wire. The fiber optic system is as described above.

Figure 49 is a flow diagram for creating a path into a vessel true lumen using laser energy, under a third alternative embodiment. Procedurally, the optical fiber system is removed from the catheter and the catheter is tracked over a standard guide wire placed in the desired sub-intimal space of the target vasculature. Because the guide wire emerges laterally from the nosecone, the tip of the catheter will track eccentrically over the guide wire.

The catheter is aligned to the vessel true lumen in accordance with Step 1, Methods 1-4 or Step 1, Method 5 above. In the case of Methods 1-4, the guide wire is retracted from the catheter, and the visualization element is advanced into the nosecone. The element is activated at the nosecone side port and the side  
5 port is rotated to face the vessel true lumen. The visualization element is removed, and the fiber optic system is re-introduced. When Method 5 is used the side port is rotated to face the vessel true lumen per fluoroscopic visualization.

The distal tip of the fiber optic system is now positioned approximately 2 centimeters proximal from the distal tip of the nosecone. At this point the  
10 application of vacuum as described above may be used to evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone. Vacuum is translated to the nosecone via the single shaft lumen.

Next, the fiber optic system is advanced distally into the nosecone. The fiber optic system is brought into contact with the sub-intimal tissue at the  
15 nosecone side port at an angle determined by the exit ramp of the nosecone. The angle of the internal ramp may vary from approximately 30 degrees to 80 degrees, but is not necessarily limited to these angles.

Next, laser energy is delivered to the optical fiber system as required to ablate the sub-intimal tissue separating the dissection plane from the vessel true  
20 lumen. If the fiber optic system does not incorporate the guide wire, the fiber optic system is removed and a standard guide wire is introduced for advancement through the pathway in the sub-intimal tissue and into the true lumen. In the case of the fiber optic system including a guide wire lumen, a guide wire is simply advanced through the lumen, through the pathway in the sub-intimal tissue, and  
25 into the true lumen.

The vacuum is now released, the guide wire position is maintained, and the catheter is retracted from the vasculature, leaving the distal portion of the guide wire in place in the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer  
5 shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; and fiber optic system outside diameter is approximately 0.010 to 0.030 inches.

#### Embodiment 5 (Method 5 under Step 3) (Figure 12)

This embodiment includes a catheter having a simple nosecone or molded  
10 distal termination, a single lumen catheter shaft, an internal slidably disposed tube, and a fiber optic system which is slidably disposed in the internal tube. The distal end of the push tube or member is slidably disposed within the catheter shaft and nosecone. Upon actuation of the push member in a distal direction, a percentage  
15 of the proximal section of the nosecone side port becomes covered. The percentage of coverage is controlled by a distal stop within the nosecone that limits the distal translation of the internal sliding member or tube. When the internal sliding member is a tube, it becomes the lumen for both the fiber optic system and guide wire.

As the internal sliding member advances into the distal position, it reduces  
20 the effective length of the nosecone side port and forces the fiber optic system or guide wire to exit the side port at a more acute angle, or an angle that is more normal to the axis of the catheter, and more normal to the sub-intimal tissue plane. This angle is governed by the proximal contact point which is against the internal sliding member, and the distal contact point which is against the exit ramp of the  
25 nosecone. The greater acute angle is designed to allow the laser energy to ablate a more direct and efficient pathway across the sub-intimal tissue.

Figures 50A and 50B show a flow diagram for creating a path into a vessel true lumen using laser energy, under a fourth alternative embodiment.

Procedurally, the sliding tube may remain in the catheter at all times. In preparation of the catheter and during delivery of the catheter to the vascular site, the distal end of the tube is retracted just proximal to the nosecone.

The optical fiber system is removed from the catheter and the catheter is tracked over a standard guide wire placed in the desired sub-intimal space of the target vasculature. Because the guide wire emerges laterally from the nosecone, the tip of the catheter will track eccentrically over the guide wire.

Next, the catheter is aligned to the vessel true lumen. This is accomplished per Step 1, Methods 1-4 or Step 1, Method 5. In the case of Methods 1-4, the guide wire is retracted from the catheter, and the visualization element is advanced into the nosecone. The element is activated at the nosecone side port and the side port is rotated to face the vessel true lumen. The visualization element is removed. In the case of Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization. The guide wire is then removed.

The internal sliding tube is then advanced until the distal end reaches the internal stop within the nosecone. The fiber optic system is introduced into the internal sliding tube and advanced to within 2 centimeters of the distal tip of the nosecone. At this point the application of vacuum can be used to evacuate fluid from the sub-intimal plane and lock the sub-intimal tissue onto the surface of the nosecone. Vacuum may be translated to the nosecone via the inner sliding tube, or via the annular lumen between the catheter shaft and the inner sliding tube. The fiber optic system is subsequently advanced distally into the nosecone, engaging the internal ramp, until a distal end reaches the exit of the side port and is in approximate contact with the sub-intimal tissue.



Laser energy is delivered to the optical fiber system as required to ablate the sub-intimal tissue separating the dissection plane from the vessel true lumen. If the fiber optic system without a guide wire is used, the fiber optic system is removed and a standard guide wire is introduced for advancement through the pathway in the sub-intimal tissue and into the true lumen. When the fiber optic system includes a guide wire lumen, a guide wire is advanced through the lumen, through the pathway in the sub-intimal tissue, and into the true lumen. The vacuum is then released, the guide wire position is maintained, and the catheter is retracted from the vasculature, leaving the distal portion of the guide wire in place in the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; internal slide tube outside diameter is approximately 0.020 to 0.030 inches; and fiber optic system outside diameter is approximately 0.010 to 0.030 inches.

#### Embodiment 6 (Method 5 under Step 3) (Figure 13)

This embodiment includes a single lumen catheter shaft, terminated in a "J" tip, used with a fiber optic system, and a standard guide wire. The "J" tip configuration can be torqued into position within the sub-intimal plane and directed towards the sub-intimal tissue. This allows the fiber optic system to be directed at an angle normal to the sub-intimal tissue, to ablate the minimum tissue necessary to establish a pathway across the sub-intimal tissue and into the vessel true lumen.

The "J" termination of the catheter can be fluoroscopically visible in an embodiment to facilitate the positioning process in the sub-intimal plane. This type of termination is fabricated or molded from fluoroscopically visible materials such

as platinum coils or gold coated stainless steel coils laminated with a variety of polymers, e.g., nylons, HDPE, and Pebax.

The "J" tip is designed to straighten in order to track over a guide wire to the vascular site, yet re-form to the "J" shape when positioned in the vasculature and the guide wire is retracted. One embodiment incorporates a visualization window just proximal to the "J" tip that is used in conjunction with an on-board visualization technique, per Step 1, Methods 1-4. For the use of IVUS, for example, this window can be fabricated from HDPE.

**Figures 51A and 51B** are a flow diagram for creating a path into a vessel true lumen using laser energy, under a fifth alternative embodiment. Procedurally, a guide wire is positioned in the sub-intimal space. The fiber optic system is removed from the catheter. The distal end of the catheter is loaded onto the guide wire. In this process the "J" tip is straightened as it tracks over the wire to the sub-intimal site. Once the terminal end of the catheter has reached the desired location, the guide wire is retracted, allowing the "J" tip to re-form.

The catheter is now aligned to the vessel true lumen per Step 1, Methods 1-4, or Step 1, Method 5. In the case of Methods 1-4, the guide wire is retracted from the catheter, and the visualization element is advanced to the visualization window. The visualization element is activated at the window and the "J" tip is rotated to face the vessel true lumen. The visualization element is then removed. In the case of Method 5, the side port is rotated to face the vessel true lumen per fluoroscopic visualization. The guide wire is then removed.

The fiber optic system is then loaded into the catheter and advanced to within 1 centimeter of the distal end of the catheter tip. Vacuum is applied through the catheter lumen per Step 2, method 1 to evacuate the sub-intimal plane. The fiber optic system is then advanced until its distal end is coincident with the distal

end of the catheter "J" tip, and is brought into approximate contact with the sub-intimal tissue.

Laser energy is delivered to the optical fiber system as required to ablate the sub-intimal tissue that separates the dissection plane from the vessel true lumen. If the fiber optic system does not include the guide wire, the fiber optic system is removed and a standard guide wire is introduced to be advanced through the pathway in the sub-intimal tissue and into the true lumen. When the fiber optic system includes a guide wire lumen, a guide wire is advanced through the lumen, through the pathway in the sub-intimal tissue and into the true lumen.

The vacuum is then released, the guide wire position is maintained, and the catheter is retracted from the vasculature, leaving the distal portion of the guide wire in place in the vessel true lumen.

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; and fiber optic system outside diameter is approximately 0.010 to 0.030 inches.

#### Embodiment 7 (Method 5 under Step 3) (Figures 14A and 14B)

This embodiment includes a catheter having a nosecone with a side port and end port connected by a slot, an internal push-ramp that is actuated by an internal push tube or member, and a fiber optic system slidably disposed in the catheter. This catheter includes a distal nosecone having a side port and a separate end port, a single lumen catheter shaft, and an internal push-tube or member that is attached to an internal hinged ramp which actuates within the nosecone side port. This catheter system may be used with any of the visualization techniques of Step 1.

The push ramp may be constructed of a flexible metal such as Nitinol or spring steel, or a polymer such as nylon or PEEK; any of which can be fabricated with or without appropriate detents to allow for bending, as required. The distal end of the push ramp is connected or hinged about the internal distal termination  
5 of the catheter shaft, opposite the side port. When the push tube or member is fully retracted proximally, the push ramp assumes a linear configuration, lying essentially flat against the inside wall of the catheter, opposite the nosecone side port.

When the push tube or member is advanced distally, the push ramp forms  
10 an incline that leads from the proximal end of the ramp, opposite the nosecone side port, to the distal end of the side port. This ramp re-directs the fiber optic system or guide wire within the nosecone side port as either is advanced distally through the nosecone side port.

Figure 52 is a flow diagram for establishing a path into a vessel true lumen  
15 using laser energy, under a sixth alternative embodiment. Procedurally, a guide wire is positioned in the sub-intimal space. The push tube or member is retracted proximally, and the distal end of the catheter is loaded onto the guide wire. The catheter is tracked over the guide wire to the desired vascular location.

The catheter is then aligned to the vessel true lumen per Step 1, Methods  
20 1-4, or Step 1, Method 5. When Methods 1-4 are used, the guide wire is retracted from the catheter, and the visualization element is advanced to the visualization window. The visualization element is activated at the nosecone side port and the side port is rotated to face the vessel true lumen. The visualization element is subsequently removed.

25 When Method 5 is used, the side port is rotated to face the vessel true lumen per fluoroscopic visualization. The guide wire is then removed.

The fiber optic system is loaded into the catheter and the distal end advanced until it is approximately 2 centimeters proximal to the nosecone or proximal to the push ramp. Vacuum is then applied through the catheter lumen per Step 2, method 1 to evacuate the sub-intimal plane. Vacuum may be applied  
5 via within the push tube or via the annular space between the push-tube and the catheter shaft. The optical fiber system is then advanced following the ramp to the nosecone cone port until it is brought into contact with the sub-intimal tissue.

Next, laser energy is delivered to the optical fiber system as required to ablate the sub-intimal tissue separating the dissection plane from the vessel true  
10 lumen. If the fiber optic system does not incorporate the guide wire, the fiber optic system is removed and a standard guide wire is introduced for advancement through the pathway in the sub-intimal tissue and into the true lumen. When the fiber optic system incorporates a guide wire lumen, a guide wire is simply advanced through the lumen, through the pathway in the sub-intimal tissue and  
15 into the true lumen.

The push tube or member is now retracted proximally, collapsing the internal push ramp to a flat configuration, and allowing the guide wire to fall through the slot and into the nosecone distal exit port. The vacuum is then released, the guide wire position is maintained, and the catheter is retracted from  
20 the vasculature. As the catheter is retracted, the guide wire will gradually translate through the slot to the nosecone distal end port, allowing the catheter to be retracted over the guide wire in a co-linear fashion, and leaving the distal portion of the guide wire in place in the vessel true lumen.

Typical dimensions of the catheter components are as follows: outer  
25 shaft/nosecone outside diameter is approximately 0.030 to 0.050 inches; internal

push tube outside diameter is approximately 0.020 to 0.030 inches; and the fiber optic system outside diameter is approximately 0.010 to 0.030 inches.

Embodiment 8 (Method 5 under Step 3) (Figure 15)

5           This embodiment includes a dual lumen catheter shaft terminating in a single lumen "J" tip, used in conjunction with an optical fiber system, an optional visualization system, and guide wires. This dual lumen catheter shaft has the same "J" type single lumen distal termination as described in Method 3, Embodiment 7, with the exception that it transitions proximally to a dual lumen  
10 catheter shaft. Fabrication and materials for the "J" tip are as described in Method 3, Embodiment 7. The dual lumen catheter shaft may be fabricated using materials and methods known in the art, and as cited in other embodiments in this description. This embodiment allows one of the slidably disposed elements included within either lumen to be advanced individually into the "J" tip single  
15 lumen.

The procedure for using this dual lumen catheter is now described with further reference to **Figure 40**. The working element described in this figure is the fiber optic system discussed below.

In a first scenario the first lumen of the dual lumen includes a standard  
20 guide wire and the second lumen includes a fiber optic system. Procedurally, the fiber optic system is loaded into its lumen and advanced until it is just proximal of the entrance to the distal single lumen ("J" tip). The catheter is loaded onto a guide wire that has been advanced into a sub-intimal plane. The catheter is then advanced over the guide wire until it reaches the desired vascular location, and the  
25 guide wire withdrawn just proximal to the entrance to the single lumen. Next, the

catheter is aligned to the vessel true lumen. This embodiment makes use of Step 1, Method 5 to align the catheter.

The sub-intimal plane is now evacuated per Step 2, Method 2, and the optical fiber system is advanced into the "J" tip so that its distal end is

5 approximately coincident with the distal tip of the catheter, or pointed towards the sub-intimal tissue. The procedure used dictates whether the tip of the fiber optic system is or is not in contact with the sub-intimal tissue. Laser energy is then delivered as required through the fiber optic system to ablate the sub-intimal tissue and create a pathway into the vessel true lumen. Once the pathway is

10 established, the distal end of the fiber optic system is retracted from the single lumen, and into the dual lumen. The standard guide wire is then advanced into the distal single lumen, and out of the "J" tip, through the pathway produced in the sub-intimal tissue, and into the vessel true lumen. Lastly, the guide wire is held in place while the catheter is retracted proximally and removed from the vasculature.

15 In a second scenario, the first lumen of the dual lumen initially contains a visualization element, per Step 1, Method 1-4, and is exchanged for the fiber optic system. The second lumen contains a standard guide wire. Procedurally, one lumen of the dual lumen shaft is loaded with the visualization element and advanced proximal to the entrance to the distal single lumen. The distal end of the

20 catheter is then loaded onto a standard guide wire that has been advanced into a sub-intimal plane. The catheter is advanced to the desired vascular location, and the guide wire is removed from the catheter and replaced with the fiber optic system.

The visualization element is then advanced into the distal single lumen

25 within the area of the window for viewing, and the "J" tip is aligned with the vessel true lumen. The visualization element is then withdrawn from the catheter, and a

standard guide wire is advance into this lumen to just proximal of the single lumen. Next, the sub-intimal plane is evacuated according to Step 2, Method 2. This may be accomplished through the lumen that houses either the guide wire, the optical fiber system, or both. Next, the optical fiber system is advanced into the "J" tip so  
5 that its distal end is approximately coincident with the distal tip of the catheter, or pointed towards the sub-intimal tissue. The tip of the fiber optic system may or may not be required to be in contact with the sub-intimal tissue, as directed by the particular procedure in use.

Laser energy is then delivered as required through the fiber optic system to  
10 ablate the sub-intimal tissue and create a pathway into the vessel true lumen.

Once the pathway has been established, the distal end of the fiber optic system is retracted from the single lumen, and into the dual lumen. The standard guide wire may then be advanced into the distal single lumen, and out of the "J" tip, through the pathway produced in the sub-intimal tissue, and into the vessel true lumen.

15 Finally, the guide wire is held in place while the catheter is retracted proximally and removed from the vasculature.

Typical dimensions of the catheter components are as follows: single lumen shaft outside diameter is approximately 0.030 to 0.050 inches; the dual lumen shaft outside diameter is approximately 0.030 to 0.050 inches (each lumen); and  
20 the fiber optic system outside diameter is approximately 0.010 to 0.030 inches.

Method 6 under Step 3: Rotational IVUS Element (Figures 9, 11-15, and  
19)

This method describes re-entry using a rotational IVUS element including a  
25 specialized distal boring tip.



**Embodiment 1 (Method 6 under Step 3) (Figure 9)**

This embodiment includes a nosecone or molded catheter termination attached to the distal end of the catheter, an internal slidably disposed actuating cannula, and an IVUS system with a specialized tip. Dimensions of the catheter components are as follows: the single lumen shaft outside diameter is approximately 0.055 inches; the cannula element outside diameter is approximately 0.040 inches; and the IVUS system outside diameter is approximately 0.030 inches.

Procedurally, with further reference to **Figure 34**, this embodiment is similar to that described for Method 3, Embodiment 3, with the exception that the IVUS system with the specialized tip provides visualization as well as the subsequent re-entry mechanism.

**Embodiment 2 (Method 6 under Step 3) (Figure 11)**

This embodiment includes a catheter having a nosecone or molded distal termination, a single lumen catheter shaft and an IVUS system with a specialized tip. Dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.045 inches; and the IVUS system outside diameter is approximately 0.030 inches.

Procedurally, with further reference to **Figure 36**, this embodiment is similar to that described for Method 3, Embodiment 5, with the exception that the IVUS system with a specialized tip provides visualization as well as the subsequent re-entry mechanism.

**Embodiment 3 (Method 6 under Step 3) (Figure 12)**

This embodiment includes a catheter having a nosecone or molded distal termination, a single lumen catheter shaft, an internal slidably disposed tube, the distal end of which translates within the catheter nosecone, and an IVUS system  
5 with a specialized tip. Dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 inches; the internal slide tube outside diameter is approximately 0.040 inches; and the IVUS system outside diameter is approximately 0.030 inches.

Procedurally, with further reference to **Figure 37**, this embodiment is similar  
10 to that described for Method 3, Embodiment 6, with the exception that the IVUS system with a specialized tip provides visualization as well as the subsequent re-entry mechanism.

**Embodiment 4 (Method 6 under Step 3) (Figure 13)**

15 The embodiment includes a single lumen catheter shaft, terminated in a "J" tip, and an IVUS system with a specialized tip. Dimensions of the catheter components are as follows: the single lumen shaft outside diameter is approximately 0.040 inches; and the IVUS system outside diameter is approximately 0.030 inches.

20 Procedurally, with further reference to **Figure 38**, this embodiment is similar to that described for Method 3, Embodiment 7, with the exception that the IVUS system with a specialized tip provides visualization as well as the subsequent re-entry mechanism.

**Embodiment 5 (Method 6 under Step 3) (Figure 14)**

This embodiment includes a catheter having a nosecone or molded distal termination, an internal push-ramp which is actuated by an internal push tube or member, a single lumen catheter shaft, and an IVUS system with specialized tip.

5 Dimensions of the catheter components are as follows: outer shaft/nosecone outside diameter is approximately 0.050 inches; the internal push tube outside diameter is approximately 0.040 inches; and the IVUS system outside diameter is approximately 0.030 inches.

Procedurally, with further reference to **Figure 39**, this embodiment is similar to that described for Method 3, Embodiment 8, with the exception that the IVUS system with a specialized tip provides visualization as well as the subsequent re-entry mechanism.

10

**Embodiment 6 (Method 6 under Step 3) (Figure 15)**

15 **Figure 15** is an embodiment that includes a dual lumen catheter shaft, terminated in a single lumen "J" tip, and an IVUS system with specialized tip. Dimensions of the catheter components are as follows: the single lumen shaft outside diameter is approximately 0.040 inches; the dual lumen shaft outside diameter is approximately 0.040 inches (per lumen); and the IVUS system outside diameter is approximately 0.030 inches.

20

Procedurally, with further reference to **Figure 40**, this embodiment is similar to that described under Method 3, Embodiment 9, with the exception that the IVUS system with a specialized tip provides visualization as well as the subsequent re-entry mechanism.

25 In general, alternatives and alternative embodiments described herein are substantially similar to previously described embodiments, and common elements

and acts or steps are identified by the same reference numbers. Only significant differences in construction or operation are described in detail. The elements and acts of the various embodiments described above can be combined to provide further embodiments.

5 All of the above references and U.S. patents and applications are incorporated herein by reference. Aspects of the invention can be modified, if necessary, to employ the systems, functions and concepts of the various patents and applications described above to provide yet further embodiments of the invention.

10 Unless the context clearly requires otherwise, throughout the description and the claims, the words "comprise," "comprising," and the like are to be construed in an inclusive sense as opposed to an exclusive or exhaustive sense; that is to say, in a sense of "including, but not limited to." Words using the singular or plural number also include the plural or singular number respectively.

15 Additionally, the words "herein," "hereunder," and words of similar import, when used in this application, shall refer to this application as a whole and not to any particular portions of this application.

The above description of illustrated embodiments of the invention is not intended to be exhaustive or to limit the invention to the precise form disclosed.

20 While specific embodiments of, and examples for, the invention are described herein for illustrative purposes, various equivalent modifications are possible within the scope of the invention, as those skilled in the relevant art will recognize. The teachings of the invention provided herein can be applied to other catheter systems, not only for the catheter system described above.

25 In general, in the following claims, the terms used should not be construed to limit the invention to the specific embodiments disclosed in the specification and

the claims, but should be construed to include all medical systems that operate under the claims. Accordingly, the invention is not limited by the disclosure, but instead the scope of the invention is to be determined entirely by the claims.

While certain aspects of the invention are presented below in certain claim  
5 forms, the inventors contemplate the various aspects of the invention in any number of claim forms. Accordingly, the inventors reserve the right to add additional claims after filing the application to pursue such additional claim forms for other aspects of the invention.

**CLAIMS**

What we claim is:

- 1     1.     A method for accessing a true lumen of a blood vessel from a sub-intimal  
2     plane of the vessel, comprising:  
3             identifying a site to enter the true lumen from a position in the sub-intimal  
4     plane distal to a chronic total occlusion (CTO);  
5             determining an orientation of the true lumen with respect to the sub-intimal  
6     plane at the selected site;  
7             physically securing tissue of the sub-intimal plane at the selected site; and  
8             establishing a path from the sub-intimal plane into the vessel true lumen.
- 1     2.     A method for crossing a chronic total occlusion (CTO) in vasculature,  
2     comprising:  
3             forming a track from a true lumen into a sub-intimal space of a blood vessel,  
4     wherein the track extends from a position proximal to the CTO in the true lumen to  
5     a position distal to the CTO in the sub-intimal space;  
6             determining an orientation of the true lumen with respect to the sub-intimal  
7     plane at an identified re-entry site from a position in the sub-intimal plane, wherein  
8     the re-entry site is distal to the CTO;  
9             physically securing tissue of the sub-intimal plane at the selected site; and  
10            selectively forming a path from the sub-intimal plane back into the true  
11    lumen.

1     3.     A catheter system for accessing a true lumen of a blood vessel from a sub-  
2     intimal plane of the vessel, comprising:  
3             at least one visualization element for determining an orientation of the true  
4     lumen with respect to the sub-intimal plane at an identified entry site from a  
5     position in the sub-intimal plane distal to a chronic total occlusion (CTO);  
6             at least one system for physically securing tissue of the sub-intimal plane at  
7     the entry-site to the catheter system; and  
8             at least one re-entry device for establishing and maintaining a path from the  
9     sub-intimal plane into the vessel true lumen.

1     4.     A catheter system for crossing chronic total occlusions (CTOs) in  
2     vasculature, comprising:  
3             means for forming a track from a true lumen into a sub-intimal space of a  
4     blood vessel, wherein the track extends from a position proximal to the CTO in the  
5     true lumen to a position distal to the CTO in the sub-intimal space;  
6             means for determining an orientation of the true lumen with respect to the  
7     sub-intimal plane at an identified re-entry site, wherein the re-entry site is distal to  
8     the CTO;  
9             means for physically securing tissue of the sub-intimal plane at the selected  
10     site; and  
11             means for selectively forming a path from the sub-intimal plane back into  
12     the true lumen.




FIGURE 1






FIGURE 2




FIGURE 3




FIGURE 4




FIGURE 5




FIGURE 6




FIGURE 7A




FIGURE 7B




FIGURE 8






FIGURE 9A




Figure 9b




FIGURE 9C




FIGURE 9D




FIGURE 9E




FIGURE 9F






FIGURE 10A








FIGURE 10D




FIGURE 10E




FIGURE 11




FIGURE 12A




FIGURE 12B




FIGURE 13




FIGURE 14A




FIGURE 14B




FIGURE 15






FIGURE 16




FIGURE 17A




FIGURE 17B




FIGURE 18




FIGURE 19




FIGURE 20




Figure 21




Figure 22






FIGURE 23




Figure 24




Figure 25




Figure 26.




Figure 27




Figure 28




FIGURE 29




FIGURE 30






FIGURE 31




FIGURE 32




FIGURE 33




FIGURE 34




FIGURE 35




FIGURE 36




FIGURE 37




FIGURE 38






FIGURE 39




FIGURE 40




FIGURE 41




FIGURE 42




FIGURE 43




FIGURE 44




FIGURE 45




FIGURE 46






FIGURE 47




FIGURE 48




FIGURE 49




FIGURE 50A




FIGURE 50B




FIGURE 51A




FIGURE 51B




FIGURE 52



(19) World Intellectual Property Organization  
International Bureau



(43) International Publication Date  
13 June 2002 (13.06.2002)


PCT

(10) International Publication Number  
**WO 02/045598 A3**

- (51) International Patent Classification<sup>7</sup>: **A61B 17/22** (71) Applicant (for all designated States except US): **LUMEND, INC.** [US/US]; 400 Chesapeake Drive, Redwood City, CA 94063 (US).
- (21) International Application Number: **PCT/US01/46871**
- (22) International Filing Date: 5 December 2001 (05.12.2001) (72) Inventors; and
- (25) Filing Language: English (75) Inventors/Applicants (for US only): **SPARKS, Kurt, D.** [US/US]; 1618 Sand Hill Road, #406, Palo Alto, CA 94304 (US). **EMERY, Jeffrey, L.** [US/US]; 2104 Meadowview Place, San Mateo, CA 94401 (US). **SEYBOLD, Brent, D.** [US/US]; 2435 Armstrong Avenue, Santa Clara, CA 95050 (US). **KUPIECKI, David, J.** [US/US]; 3276 Market Street, San Francisco, CA 94114 (US). **PINSON, C., Danielle** [US/US]; 711 Eunice Avenue, Mountain View, California 94040 (US). **MADSEN, Allen, W.** [US/US]; 2878 Kilo Avenue, San Jose, CA 95124 (US). **KELEHER, Michael, D.** [US/US]; 34330 Dobson Way, Fremont, CA 94555 (US). **SALINAS, Sergio** [US/US]; 642 Scott Avenue, Redwood City, CA 94063 (US). **CLARK, Benjamin, J.** [US/US]; 3814 East Lake Way, Redwood City, CA 94062 (US). **SELMON, Matthew, R.** [US/US]; 99 Walnut Avenue, Atherton, CA 94027 (US).
- (26) Publication Language: English
- (30) Priority Data:
- |            |                               |    |
|------------|-------------------------------|----|
| 60/251,756 | 5 December 2000 (05.12.2000)  | US |
| 60/255,729 | 14 December 2000 (14.12.2000) | US |
| 60/263,350 | 22 January 2001 (22.01.2001)  | US |
| 60/263,397 | 22 January 2001 (22.01.2001)  | US |
| 60/263,579 | 22 January 2001 (22.01.2001)  | US |
| 60/263,580 | 22 January 2001 (22.01.2001)  | US |
| 60/263,589 | 22 January 2001 (22.01.2001)  | US |
| 60/268,263 | 12 February 2001 (12.02.2001) | US |
| 60/301,537 | 27 June 2001 (27.06.2001)     | US |
| 60/329,936 | 17 October 2001 (17.10.2001)  | US |

[Continued on next page]

(54) Title: CATHETER SYSTEM FOR VASCULAR RE-ENTRY FROM A SUB-INTIMAL SPACE



(57) Abstract: A catheter system and corresponding methods are provided for accessing a blood vessel true lumen from a sub-intimal plane of the vessel. The catheter system includes visualization elements for determining the orientation of the true lumen with respect to the sub-intimal plane at an identified entry site from a position in the sub-intimal plane. The entry site is distal to a chronic total occlusion (CTO). The catheter system also includes a system for physically securing tissue of the sub-intimal plane at the entry site to the catheter system. The attaching system reduces or eliminates catheter float within the sub-intimal space. The catheter system further includes re-entry devices to establish and maintain a path from the sub-intimal plane back into the vessel true lumen.



(74) **Agents:** GREGORY, Richard, L. et al.; Perkins Coie LLP, P.O. Box 2168, Menlo Park, CA 94026 (US).

European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

(81) **Designated States (national):** AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZM, ZW.

**Published:**

— with international search report

(88) **Date of publication of the international search report:**

16 January 2003

(84) **Designated States (regional):** ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM),

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

Internal Application No.

PCT/US 01/46871

### A. CLASSIFICATION OF SUBJECT MATTER

IPC 7 A61B17/22

According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 7 A61B

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

EPO-Internal

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 99 35980 A (LUMEND INC ;MILO CHARLES F (US); VETTER JAMES W (US); HINOHARA TOM) 22 July 1999 (1999-07-22) page 15, line 1 - line 5; figures 5,6 page 17, line 29 -page 18, line 17; figure 4	3
Y	---	4
Y	WO 99 35979 A (CO FRED H ;FRENCH RONALD G (US); HILL RICHARD E (US); LUMEND INC ( ) 22 July 1999 (1999-07-22) abstract; figures 3-4D	4
A	---	3
A	US 6 068 638 A (MAKOWER JOSHUA) 30 May 2000 (2000-05-30) column 4, line 6 -column 5, line 8; figure 2	3,4
	---	

☒ Further documents are listed in the continuation of box C.

☒ Patent family members are listed in annex.

\* Special categories of cited documents :

- \*A\* document defining the general state of the art which is not considered to be of particular relevance
- \*E\* earlier document but published on or after the International filing date
- \*L\* document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- \*O\* document referring to an oral disclosure, use, exhibition or other means
- \*P\* document published prior to the International filing date but later than the priority date claimed

- \*T later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- \*X document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- \*Y document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.
- \*A document member of the same patent family

Date of the actual completion of the international search

13 September 2002

Date of mailing of the international search report

23/09/2002

Name and mailing address of the ISA  
European Patent Office, P.B. 5818 Patentlaan 2  
NL - 2280 HV Rijswijk  
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,  
Fax: (+31-70) 340-3016

Authorized officer \_\_\_\_\_

Ducreau, F

# INTERNATIONAL SEARCH REPORT

Interj 1st Application No  
PCT/US 01/46871

C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT		
Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	US 5 951 541 A (GERTNER KEVIN ET AL) 14 September 1999 (1999-09-14) column 8, line 8 - line 62; figures 2-4 -----	3,4

# INTERNATIONAL SEARCH REPORT

ational application No.  
PCT/US 01/46871

## Box I Observations where certain claims were found unsearchable (Continuation of item 1 of first sheet)

This International Search Report has not been established in respect of certain claims under Article 17(2)(a) for the following reasons:

1. ☒ Claims Nos.: 1,2  
because they relate to subject matter not required to be searched by this Authority, namely:  
Rule 39.1(iv) PCT - Method for treatment of the human or animal body by surgery
2. ☐ Claims Nos.:  
because they relate to parts of the International Application that do not comply with the prescribed requirements to such an extent that no meaningful International Search can be carried out, specifically:
3. ☐ Claims Nos.:  
because they are dependent claims and are not drafted in accordance with the second and third sentences of Rule 6.4(a).

## Box II Observations where unity of invention is lacking (Continuation of item 2 of first sheet)

This International Searching Authority found multiple inventions in this International application, as follows:

1. ☐ As all required additional search fees were timely paid by the applicant, this International Search Report covers all searchable claims.
2. ☐ As all searchable claims could be searched without effort justifying an additional fee, this Authority did not invite payment of any additional fee.
3. ☐ As only some of the required additional search fees were timely paid by the applicant, this International Search Report covers only those claims for which fees were paid, specifically claims Nos.:
4. ☐ No required additional search fees were timely paid by the applicant. Consequently, this International Search Report is restricted to the invention first mentioned in the claims; it is covered by claims Nos.:

Remark on Protest

- ☐ The additional search fees were accompanied by the applicant's protest.
- ☐ No protest accompanied the payment of additional search fees.

# INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/US 01/46871

Patent document cited in search report		Publication date	Patent family member(s)	Publication date
WO 9935980	A	22-07-1999	AU 2004999 A	02-08-1999
			WO 9935980 A1	22-07-1999
WO 9935979	A	22-07-1999	US 6231546 B1	15-05-2001
			AU 745116 B2	14-03-2002
			AU 1831599 A	02-08-1999
			BR 9814588 A	24-10-2000
			CA 2318470 A1	22-07-1999
			EP 1047344 A1	02-11-2000
			JP 2002508987 T	26-03-2002
			WO 9935979 A1	22-07-1999
			US 6221049 B1	24-04-2001
			US 6235000 B1	22-05-2001
US 6068638	A	30-05-2000	US 5830222 A	03-11-1998
			US 6231587 B1	15-05-2001
			AU 1275997 A	11-05-1998
			AU 733332 B2	10-05-2001
			AU 1847297 A	22-08-1997
			AU 723785 B2	07-09-2000
			AU 1847397 A	22-08-1997
			CA 2244079 A1	07-08-1997
			CA 2244080 A1	07-08-1997
			CN 1218414 A	02-06-1999
			CN 1216929 A	19-05-1999
			EP 0879068 A1	25-11-1998
			EP 0932426 A1	04-08-1999
			JP 2000504594 T	18-04-2000
			JP 2000505316 T	09-05-2000
			US 6283951 B1	04-09-2001
			WO 9816161 A1	23-04-1998
			WO 9727893 A1	07-08-1997
			WO 9727898 A1	07-08-1997
			US 6159225 A	12-12-2000
			US 6283983 B1	04-09-2001
			US 6375615 B1	23-04-2002
			US 2002002349 A1	03-01-2002
			US 6302875 B1	16-10-2001
			US 2002029079 A1	07-03-2002
			US 2001047165 A1	29-11-2001
			US 6432127 B1	13-08-2002
			US 2002062146 A1	23-05-2002
			AU 726713 B2	16-11-2000
			AU 7431696 A	30-04-1997
			AU 729466 B2	01-02-2001
			AU 7595396 A	30-04-1997
			CA 2234361 A1	17-04-1997
			CA 2234389 A1	17-04-1997
			EP 1166721 A2	02-01-2002
			EP 0954248 A1	10-11-1999
			EP 0955933 A1	17-11-1999
			EP 0910298 A1	28-04-1999
			JP 11513577 T	24-11-1999
			JP 11514269 T	07-12-1999
			WO 9713471 A1	17-04-1997
			WO 9713463 A1	17-04-1997
			US 6190353 B1	20-02-2001

# INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/US 01/46871

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
US 5951541	A	14-09-1999	NONE

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

☒ **BLACK BORDERS**

☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**

☐ **FADED TEXT OR DRAWING**

☒ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**

☐ **SKEWED/SLANTED IMAGES**

☒ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**

☐ **GRAY SCALE DOCUMENTS**

☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**

☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**

☐ **OTHER:** \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**